

Corning® Centrifuge Tubes

Designed for Optimum Performance

CORNING

Corning offers a full line of disposable, premium quality conical centrifuge tubes manufactured from high quality, ultra-clear resins. These tubes meet the quality and performance standards demanded by researchers.

► **High Speed 15 mL and 50 mL tubes**

- Rated up to 17,000 RCF
- 95 kPA (14 psi) pressure-tested (plug seal cap)
- RNase-/DNase-free
- Nonpyrogenic
- Sterile

► **Large Volume 250 mL and 500 mL tubes**

- Rated up to 6,000 RCF
- Best choice for cell harvesting and purification
- Nonpyrogenic
- Sterile

► **Simplify your workflow with:**

- Color bands visually identify samples
- Hazardous materials pictograms that comply with OSHA and GHS standards



Ordering Information



15 mL and 50 mL centrifuge tubes with CentriStar cap



Easy ID color strip tubes



50 mL centrifuge tube with hazardous materials label

Corning® 15 mL Centrifuge Tubes

Cat. No.	Material	Cap Style	Max. RCF	Feature	Qty/Pk	Qty/Cs
430053	PET	Plug seal	3,600	Conical bottom	25/sleeve	500
430055	PET	Plug seal	3,600	Conical bottom	50/rack	500
430052	PP	Plug seal	12,000	Conical bottom	50/rack	500
430766	PP	Plug seal	12,000	Conical bottom	25/sleeve	500
430790	PP	CentriStar™	12,500	Conical bottom	50/rack	500
430791	PP	CentriStar	12,500	Conical bottom	25/sleeve	500

Corning 50 mL Centrifuge Tubes

430304	PET	Plug seal	3,600	Conical bottom	25/rack	500
430290	PP	Plug seal	15,500	Conical bottom	25/rack	500
430291	PP	Plug seal	15,500	Conical bottom	25/sleeve	500
430828	PP	CentriStar	17,000	Conical bottom	25/rack	500
430829	PP	CentriStar	17,000	Conical bottom	25/sleeve	500
4558	PP	CentriStar	17,000	Conical bottom	25/Universal rack*	300
New 431526	PP	Plug seal	15,500	Easy ID color strip - green, conical bottom	25/Universal rack	300
New 431527	PP	Plug seal	15,500	Easy ID color strip - Orange, conical bottom	25/Universal rack	300
New 431528	PP	Plug seal	15,500	Easy ID color strip - blue, conical bottom	25/Universal rack	300
New 431525	PP	Plug seal	15,500	Hazardous materials label, conical bottom	25/Universal rack	300
430897	PP	Plug Seal	3,000	Self-standing bottom	25/sleeve	500
430921	PP	CentriStar	3,000	Self-standing bottom	25/sleeve	500

*New innovative universal rack can hold 50 mL and 15 mL tubes securely, allowing researchers to work with and store both size tubes in the same rack, saving bench and storage space.

250 mL and 500 mL Premium Quality Centrifuge Tubes and Support Cushions

Cat. No.	Description	Material	Cap Style	Max RCF	Qty/Pk	Qty/Cs
430776	250 mL Tube	PP	Plug	6,000	6	102
430236	250 mL Support cushion*	PEI	N/A	N/A	N/A	6
431123	500 mL Tube	PP	Plug	6,000	6	36
431124	500 mL Support cushion*	PEI	N/A	N/A	n/a	6

*Support cushions must be used with the 250 mL and 500 mL centrifuge tubes unless the rotor has appropriately shaped V-bottom holders. PP = Polypropylene, PET = Polyethylene Terephthalate, RCF = Relative Centrifugal Force (x g), PEI = Polyetherimide

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

Corning Incorporated Life Sciences

836 North St.
Building 300, Suite 3401
Tewksbury, MA 01876
t 800.492.1110
t 978.442.2200
f 978.442.2476

www.corning.com/lifesciences

Worldwide Support Offices

ASIA/PACIFIC
Australia/New Zealand
t 0402-794-347

China
t 86 21 2215 2888
f 86 21 6215 2988

India
t 91 124 4604000
f 91 124 4604099

Japan
t 81 3-3586 1996
f 81 3-3586 1291

Korea
t 82 2-796-9500
f 82 2-796-9300

Singapore
t 65 6733-6511
f 65 6861-2913

Taiwan
t 886 2-2716-0338
f 886 2-2516-7500

EUROPE

France
t 0800 916 882
f 0800 918 636

Germany
t 0800 101 1153
f 0800 101 2427

The Netherlands
t 31 20 655 79 28
f 31 20 659 76 73

United Kingdom
t 0800 376 8660
f 0800 279 1117

All Other European Countries

t 31 (0) 20 659 60 51
f 31 (0) 20 659 76 73

LATIN AMERICA Brasil

t (55-11) 3089-7419
f (55-11) 3167-0700

Mexico
t (52-81) 8158-8400
f (52-81) 8313-8589

CORNING | **FALCON** | **AXYGEN** | **GOSELIN** | **PYREX**

For a listing of trademarks, visit us at www.corning.com/lifesciences/trademarks.
All other trademarks in this document are the property of their respective owners.