

# Break the Winter Cold with these HOT SPECIALS from LPS!

## EQUIPMENT



### RotoBot® Mini Programmable Rotator

- Offers advance programming with a wide variety of rotational modes
- Cold room and incubator safe
- Accommodates up to 50mL tubes
- Includes 6x1.5mL holder (2ea.), 1x15mL holder (2ea.) and 1 x 50mL holder (2 ea.)

Cat. #	Description	List	SALE
R4040	RotoBot Mini Programmable Tube Rotator, with starter set	\$495.00	<b>\$371.25</b>

### Prism™ Air-Cooled Microcentrifuge

- Unique multi-flow air cooling technology
- Updated control panel, easy to read LCD
- Includes 24 place (15000rpm/21,200xg) rotor
- Quick acceleration and deceleration time



Cat. #	Description	List	SALE
C2500	Prism™ Microcentrifuge	\$1,816.15	<b>\$1,450.00</b>



### ISOBlock® Digital Dry Baths

- Two isolated chambers with dual temperature control
- Available QuickFlip blocks for microtubes & PCR plates
- Blocks not included
- Hinged lid improves temperature uniformity
- Ambient +5 to 105°C temperature range

Cat. #	Description	List	SALE
BSH6000	IsoBlock® Digital Dry Bath	\$659.00	<b>\$494.25</b>
BSW1500	Block, 24x1.5mL (conical) Centrifuge Tubes	\$98.00	<b>\$73.50</b>
BSWCMB	Quick-Flip Block, holds 24x1.5mL or 32x0.2mL & 20x.5mL Tubes	\$129.00	<b>\$96.75</b>

### Accuris UV Transilluminator

- 302nm midrange UV wavelength
- Dual intensity settings (50% or 100%)
- 16x19cm viewing area
- Hinged UV safety cover



Cat. #	Description	List	SALE
E3000	Accuris UV Transilluminator	\$930.00	<b>\$790.50</b>

Buy the E3000 and get the E5000-SD Smart Doc Enclosure

**FREE**  
(a \$550.00 value)



### The Edge™ Dedicated pH/ORP Tablet Meter

- The world's most innovative pH Meter
- Large easy to read LCD
- Capacitive touch for on-screen control
- Incredibly thin! Use as a wall mount or cradle mount
- Includes meter, temp probe, electrode, cradle, pH packets



Cat. #	Description	List	SALE
HI 2002	Edge™ Dedicated pH/ORP Meter	\$391.00	<b>\$332.35</b>

# SUPPLIES

## PCR PLATES


### Nest® 96 well PCR Plates



- Virgin thin walled polypropylene promotes fast thermal transfer
- Available in semi-skirted and non-skirted versions
- Available with elevated wells to minimize cross contamination
- White plates are ideal for qPCR applications

Cat. #	Description	List	SALE
402001	PCR Plate, non-skirted, Clear 200uL, 25/pk, 100/cs	\$216.00	<b>\$150.00</b>
402001-A	PCR Plate, non-skirted, Clear 200uL, 10/pk	\$21.60	<b>\$17.28</b>
402301	PCR Plate, semi-skirted, Clear 200uL, 25/pk, 100/cs	\$216.00	<b>\$150.00</b>
402301-A	PCR Plate, semi-skirted, Clear 200uL, 10/pk	\$21.60	<b>\$17.28</b>
402401	PCR Plate, semi-skirted, Clear 100uL low profile, 25/pk, 100/cs	\$216.00	<b>\$150.00</b>


### Thermaseal2 PCR Film



- Recommended range of -40° to 120°C
- Certified RNase, DNase and nucleic-acid free

Cat. #	Description	List	SALE
T329020	ThermalSeal2 PCR Film, 2mil thick, non-sterile, 100/pk	\$72.00	<b>\$57.60</b>


### 384-Well Thin Wall PCR Plates



- Designed for high throughput laboratory work
- Skirted to allow bar-coding on side for identification
- Notched corners are cut away to facilitate proper orientation
- Compatible with automated fluid handling systems or multichannel pipettors
- Each well has 40µl capacity
- Opaque plates ensure low levels of background fluorescents

Cat. #	Description	List	SALE
L234350	384-well PCR Plate, full skirt, A24, One Notch, natural, 100/cs	\$463.83	<b>\$295.00</b>

### AlumaSeal 96 PCR Plate Foil



- Recommended for 96-well plates
- Heat and cold resistant, recommended for temperatures from -80°C to +150°C
- Certified DNase-, RNase-, and nucleic-acid-free
- 38µm thick aluminum foil film fits inside the raised rim plates

Cat. #	Description	List	SALE
T329096	AlumaSeal 96 Foils non-sterile, 100/pk	\$63.24	<b>\$50.60</b>

## CENTRIFUGE TUBES

### SpectraTube® Rainbow Cap Centrifuge Tubes

- Color code your workflow or just brighten up the lab
- Certified DNase, RNase, & DNA free
- Max g-Force of 12,000 (50mL) and 15,000xg (15mL)
- Five assorted colors (blue, green, yellow, red, orange) in foam racks
- Also available in bulk packaging

Cat. #	Description	List	SALE
L232715	SpectraTube® with Assorted Caps, 15mL, Str1, 500/cs	\$139.00	<b>\$89.00</b>
L232750	SpectraTube® with Assorted Caps, 50mL, Str1, 500/cs	\$159.00	<b>\$99.00</b>



# SURE STRAIN PREMIUM CELL STRAINERS

## SureStrain™ Premium Cell Strainers



- Advanced design for faster flow through
- Superior ventilation prevents clogging and spills
- Sterile, RNase, DNase, pyrogen and cytotoxin free
- Packed in individual blister packs

Cat. #	Description	List	SALE
M854040	SureStrain™ Premium Cell Strainers, 40µm, Blue, 50/pk	\$99.00	<b>\$79.20</b>
M854070	SureStrain™ Premium Cell Strainers, 70µm, White, 50/pk	\$99.00	<b>\$79.20</b>
M854100	SureStrain™ Premium Cell Strainers, 100µm, Yellow, 50/pk	\$99.00	<b>\$79.20</b>

## SureStrain™ Reducing Adapter



- Enables the cell strainers to be used with 15mL and 5mL tubes
- Packed individually sterile, 25/pk

Cat. #	Description	List	SALE
M854015	SureStrain™ Reducing Adapter for 5mL & 15mL tubes, 25/pk	\$69.00	<b>\$55.20</b>

## UPrep® Spin Columns

- Silica spin basket in 2 mL tube
- Ideal for use in plasmid miniprep and gel extraction kits
- Use w/manufacturers PCR clean up & genomic DNA kits
- No change in protocols or yields, works w/most kits



Cat. #	Description	List	SALE
M858143	UPrep® Spin Filters with basket and tubes, 50/unit	\$44.87	<b>\$35.90</b>

# TUBE RACKS


## Adapt-A-Rack® 4 Place Combo Tube Rack



- Flexible multi tube rack holds all 5-50mL tubes
- Open sides allow for label viewing
- Upright racks can be linked together as needed
- Available in a combination of colors, 2 per pack

Cat. #	Description	List	SALE
HS120184	Adapt-A-Rack Combo, holds 4 tubes 5-50mL, Blue/Green, 2 ea	\$13.90	<b>\$10.43</b>
HS120185	Adapt-A-Rack Combo, holds 4 tubes 5-50mL, Blue/Yellow, 2 ea	\$13.90	<b>\$10.43</b>

## Snap-Together Conical Tube Rack



- Orange ABS rack holds 30 x15mL and 20x50 mL conical tubes
- Compact design for stable use in water baths
- Alpha-numeric grid for easy identification
- Measures 8 x 6.75 x 2.4"H, shipped flat, easy to assemble

Cat. #	Description	List	SALE
HS24306	Snap-Together Conical Tube Rack, 30x15 & 20x50mL, ea	\$25.72	<b>\$19.29</b>

## Pop-Up™ Conical Tube Racks



- Closes flat for compact storage, pops up for quick use
- Simple, easy action to open or close the rack
- Imprinted grid helps organize 15 and 50 mL conical tubes
- Available in dedicated or multi-tube configurations
- Made of rugged polypropylene, 2 per pack

Cat. #	Description	List	SALE
T524320	Pop-Up Rack holds 21 x 15mL & 12 x 50 mL, Pink, 2/pk	\$35.75	<b>\$26.82</b>



## Laboratory Products Sales Inc.

*A NY State Certified Woman-Owned Small Business*

1665 Buffalo Road • Rochester, N.Y. 14624 • 585-247-4720

TOLL FREE: 1-800-388-0166 • FAX: 585-247-6686 • www.LPSinc.com