

Prices and offers valid from 1/1/10 through 3/31/10

LPS SAVES YOU MORE \$\$\$ ON YOUR LIQUID HANDLING SUPPLIES

MICROCENTRIFUGE TUBES

Cat. #	Description	Regular	Sale Price
L211511-CS	1.7mL MicroCentrifuge Tubes, grad, natural, RNase/DNase free, 5,000/case	\$140.90	\$80.00



**PRICING
BY THE CASE
SAVES YOU
MORE!**

PCR TUBES INDIVIDUAL TUBES

- Manufactured of autoclavable virgin polypropylene
- Designed for optimal fit in all leading thermal cyclers
- Thin walls accelerate even heat transfer and shorten cycle times
- Independently tested and certified RNase, DNase free

Cat. #	Description	Regular	Sale Price
L243001-CS	200µL thin wall PCR tube w/dome caps, 5 bags of 1000 per case	\$237.25	\$136.52
L243101-CS	200µL PCR tube, w/attached flat caps, 5 bags of 1000 per case	\$207.40	\$142.30



**PRICING
BY THE CASE
SAVES YOU
MORE!**

PCR STRIP TUBES

- Thin-wall design provides maximum PCR efficiency
- Easily cut for use as individual tubes
- Working volume = 200µL

Cat. #	Description	Regular	Sale Price
L216380-CS	200µL 8-strip PCR tube w/separate dome caps, 960 per unit, 5 units per case	\$440.80	\$325.00
L216382-CS	200µL 8-strip PCR tube w/separate flat caps, 960 per unit, 5 units per case	\$478.00	\$325.00



SCREW CAP MICROCENTRIFUGE TUBES

- Suitable for vapor phase LN2 storage
- Temperature range of -190°C to +121°C
- Clear polypropylene for sample ID and chemical resistance
- Available sterile (w/caps) or non-sterile (w/o caps)



STERILE WITH ATTACHED CAPS, 1000/UNIT

Cat. #	Description	Regular	Sale Price
L233032	0.5mL o-ring tube, skirted	\$188.00	\$108.87
L233042	1.5mL o-ring tube, conical	\$188.00	\$108.87
L233052	1.5mL o-ring tube, skirted	\$188.00	\$108.87
L233062	2.0mL o-ring cap tube, skirted, conical	\$188.00	\$108.87
L233072	2.0mL o-ring cap tube, skirted, conical	\$188.00	\$108.87

NON-STERILE WITHOUT CAPS, 1000/UNIT

Cat. #	Description	Regular	Sale Price
L233021	0.5mL tube, conical	\$48.05	\$28.00
L233031	0.5mL tube, skirted, conical	\$48.05	\$28.00
L233033	0.5mL tube, skirted, conical, Brown	\$72.00	\$51.00
L233041	1.5mL tube, conical	\$48.05	\$28.00
L233043	1.5mL tube, conical, Brown	\$72.00	\$51.00
L233051	1.5mL tube, skirted, conical	\$48.05	\$28.00
L233053	1.5mL tube, skirted, conical, Brown	\$72.00	\$51.00
L233061	2.0mL tube, conical	\$48.05	\$28.00
L233063	2.5mL tube, conical, Brown	\$72.00	\$51.00
L233071	2.0mL tube, skirted, conical	\$48.05	\$28.00
L233073	2.0mL tube, skirted, conical, Brown	\$72.00	\$51.00

SCREW CAPS WITH O-RINGS, 1000/UNIT

Cat. #	Description	Regular	Sale Price
L233121	Screw cap with O-ring Natural, PP	\$60.47	\$44.64
L233122	Screw cap with O-ring Red, PP	\$60.47	\$44.64
L233123	Screw cap with O-ring Blue, PP	\$60.47	\$44.64
L233124	Screw cap with O-ring Green, PP	\$60.47	\$44.64
L233125	Screw cap with O-ring Yellow, PP	\$60.47	\$44.64
L233126	Screw cap with O-ring Orange, PP	\$60.47	\$44.64
L233127	Screw cap with O-ring Dark Brown, PP	\$60.47	\$44.64
L233128	Screw cap, with O-Ring, Violet, PP	\$60.47	\$44.64
L233129	Screw cap with O-ring Assorted, PP	\$60.47	\$44.64

SNAP CAP STERILE CULTURE TUBES

- Precision molded polypropylene (PP)
- Supplied with two-position ribbed polyethylene caps for aerobic or anaerobic closing
- Polypropylene tubes tolerate higher centrifugation speeds (approx. 1,000 xg) and temperatures from -190°C to 121°C
- NOTE: The individually wrapped tubes have printed graduations and writing area.

5mL (12 x 75mm)

Cat. #	Description	Regular	Sale Price
L240502	12x75mm tubes, PP, 5mL, STRL, indiv wrapped, 500/CS	\$93.91	\$77.00
L285681	12x75mm tubes, PP, 5mL, STRL, 125/bg, 1M/CS	\$102.17	\$70.00
L285861	12x75mm tubes, PP, 5mL, STRL, 25/pk, 500/CS	\$55.61	\$41.13

14mL (17 x 100mm)

Cat. #	Description	Regular	Sale Price
L240602	17x100mm tubes, PP, 14mL, STRL, indiv wrapped, 500/CS	\$116.59	\$96.98
L285991	17x100mm tubes, PP, 14mL, STRL, 25/bg, 500/CS	\$90.29	\$59.00



REAGENT RESERVOIRS

- Sharp -V- Bottom eliminates waste
- RNase/DNase/DNA Free and non-pyrogenic
- Unique stacking "lugs" for easy removal from stack
- Graduated with 4 pour off corners

Cat. #	Description	Regular	Sale Price
M255100	Reagent reservoir, 25mL Sterile, indiv wrap, 100/cs	\$63.00	\$48.51
M255200	Reagent reservoir, 25mL Sterile, 5/pack, 200/cs	\$80.00	\$61.60
M501100	Reagent reservoir, 50mL Sterile, indiv wrap, 100/cs	\$63.00	\$45.00
M505200	Reagent reservoir, 50mL Sterile, 5/pack, 200/cs	\$80.00	\$61.60
M100110	Reagent reservoir, 100mL Sterile, indiv wrap, 100/cs	\$56.70	\$45.00
M100520	Reagent reservoir, 100mL Sterile, 5/pack, 200/cs	\$80.00	\$61.60



Laboratory Products Sales Inc.

1665 Buffalo Road • Rochester, N.Y. 14624 • 585-247-4720

TOLL FREE: 1-800-388-0166 • FAX: 585-247-6686 • www.LPSinc.com