

# LAB VALUES

from  
*Labnet International*

Winter 2019


AGITATION  
CONSTANT TEMPERATURE

Special offers valid - January 1 - March 31, 2019

## Rocker 35 EZ



- Large capacity rocker
- Can easily handle up to 11 lbs of load
- Tilt easily adjustable with knob on the side
- Single and dual platform models
- **5-yr warranty**



Single Platform \$618  
Double Platform \$680

### SPECIFICATIONS

Speed range 10-120 rpm  
Timer 10 - 120 min or continuous  
Tilt range Adjustable,  $\pm 20^\circ$   
Ambient operating range  $-10^\circ$  to  $70^\circ\text{C}$ , up to 85%RH

WLV-1.19

### Ordering Information


S2035-E *Rocker 35 EZ*, single platform  
S2035-DE *Rocker 35 EZ*, double platform

- Left to right linear action
- Gentle to vigorous shaking
- Four platform configurations
- Safe for use in temperature controlled environments

### SPECIFICATIONS

Speed range 6-150 strokes per min  
Timer 20- 120 min or continuous  
Motion Linear, 30 mm stroke length  
Ambient operating range  $+4^\circ$  to  $70^\circ\text{C}$

## Pro 30 Reciprocal Shaker



### Ordering Information

S2030-RC-B	<i>Pro 30 Reciprocal Shaker</i> , without platform	<b>\$692</b>
S2031-12	Flat platform with non-slip rubber mat	\$119
S2031-12D	As above, but two stacked platforms	\$194
S2031-13	Spring loaded platform for flasks, bottle or tube racks	\$221
S2031-18	Flask clamp platform without clamps, pre-drilled	\$194

WLV-2.19

SAVE A COOL 20%

See Page 4  
for Lab Values  
on other  
product lines

# VX-200 Vortex Mixer



- CombiCup™ accepts variety of tube sizes
- Circular orbit for effective vortexing at any speed
- Full range of accessories

## SPECIFICATIONS

Speed range	0-3,400 rpm
Timer	Touch or continuous
Ambient operating range	+4° to 65°C

## Ordering Information

S0200	VX-200 Vortex Mixer with Combi Cup Head	<b>\$141</b>			
<b>Optional head attachments</b>					
S0200-20	3" flat head adapter with dimpled mat	\$28	S0200-24	6 x 50 mL tubes	\$34
S0200-21	24 x 1.5/2.0 mL, 24 x 0.5, and 32 x 0.2 mL tubes	\$28	S0200-25	12 x 1.5/2.0 mL tubes, held horizontally	\$34
S0200-22	1 microplate or 64 x 0.2 mL tubes or 8 x 0.2 mL tube strips	\$34	S0200-26	4 x 15 mL tubes, held horizontally	\$66
S0200-23	8 x 15 and 8 x 12/13 mm diameter tubes	\$34	S0200-27	2 x 50 mL tubes, held horizontally	\$66

WLV-3-19

- Single, dual and four block capacity
- Dual display -- time & temperature
- USB connectivity enables traceability of data
- Molded block chamber
- Uses same blocks - same reliable performance

## Accublock Digital Dry Bath & Dry Bath Blocks



Save 20% on blocks



**D1301**  
8.3 x 12.4 x 4.7 in  
**\$283**



**D1302**  
8.3 x 12.4 x 4.7 in  
**\$364**



**D1304**  
8.3 x 16.3 x 4.7 in  
**\$538**

## Ordering Information

- D1301 **Digital Dry Bath**, single block capacity  
 D1302 **Digital Dry Bath**, dual block capacity  
 D1304 **Digital Dry Bath**, four block capacity

WLV-4.19


# Space Saving Equipment for your Lab Bench



Mini Dry Bath  
5.5 x 4.8 x 2.4 in  
**\$218**



Mini Incubator  
11 x 11 x 3.2 in  
**\$314**



Gyro Mini 3-D Shaker  
10 x 7 x 6 in  
**\$239**

## Ordering Information

- S0500 *GyroMini 3-D shaker, with dimpled rubber mat*  
 D0100 *AccuBlock Mini Dry Bath, with 12 x 1.5 mL block*  
 I5110A *Mini Incubator*

WLV-5.19

**HEATS & COOLS**

# AccuTherm Shaking Incubator



Block B



Block C



Block D



Block H



Block J

## Ordering Information

- |            |  |                |
|------------|--|----------------|
| I-4001-HCS | <i>AccuTherm Microtube Shaking Incubator</i> | <b>\$1,598</b> |
| I-4000-B   | Block B, 54 x 0.5mL                          | <b>\$ 151</b>  |
| I-4000-C   | Block C, 96 x 0.2mL                          | <b>\$ 151</b>  |
| I-4000-D   | Block D, 24 x 15mm                           | <b>\$ 151</b>  |
| I-4000-H   | Block H, 40 x 2.0mL                          | <b>\$ 170</b>  |
| I-4000-J   | Block J, 96 Well ELISA Plate                 | <b>\$ 231</b>  |

WLV-6.19

## FastPette™ V2 & FastPette™ Pro Pipette Controllers

### FastPette V2

- Light weight handle with smooth push-button
- Ergonomically located switches
- Works with glass & plastic pipettes in the 0.5 - 100 range
- Fully charged can be used continuously for 8 hours

**\$239**  
~~\$299~~



### FastPette Pro

- Convenient ergonomically located switches
- Designed for use with glass or plastic serological pipettes from 0.5 to 100mL range
- NiMH batteries when fully charged - 8 hours continuous use
- LCD screen provides control over the battery status

**\$300**  
~~\$375~~



### Ordering Information

P2000 *Fastpette V2* Pipette Controller

P2002 *Fastpette Pro* Pipette Controller

WLV-7.19

## Labpette R™ & Combi-Syringe Tips

- Ideal for repetitive dispensing of all difficult liquids
- Convenient tip selection chart located on body of pipette
- Maintenance free



**\$274**



**\$80**

### Ordering Information

P3000 *Labpette R Repeating Pipette*, 1 µl to 5,000 µl

P3535 *Combination pack*, 20 each of 0.5, 1.25, 2.5, 5.0 and 12.5 mL tips

WLV-8.19

## High Performance Tabletop REFRIGERATED Centrifuge



- Large graphic LCD screen
- Dedicated button for fast cooling
- Maintenance free induction motor
- 7 rotor options
- 99 user programs

**\$5,559**



All rotors and accessories are discounted **20%**. Check with your dealer for pricing and part #'s

### Ordering Information

C0226R *High Performance Refrigerated Centrifuge*

WLV-9.19

## High Performance Benchtop Centrifuge

- 2 models-ambient & refrigerated
- Digital display and safety indicators
- Angled front panel and large knobs
- Optional external temperature controller



**\$3,178**  
C0336



**\$4,959**  
C0336R

All rotors and accessories are discounted **20%**. Check with your dealer for pricing and part #'s

### Ordering Information

C0336 *High Performance Benchtop Centrifuge*, ambient

C0336R *High Performance Benchtop Centrifuge*, refrigerated

## Enduro UV Transilluminators

- Designed for rugged use
- Compact size and quality offer a great value
- Single and dual wavelength models

Single wave length \$1,213  
Dual wave length \$1,281



### Ordering Information

U1001 *UV Transilluminator with 302 nm wave length*

U1002 *UV Transilluminator with 302 and 365 nm wavelengths*

WLV-11.19

## Enduro™ GDS & GDST\*

\*Gel Documentation System

- Plug & play installation <15 min
- Small footprint
- 4 click image acquisition
- Integrated system (GDST)

**\$4,800**  
(List \$6,272)

Laptop not included



### Ordering Information

GDS-1302 *GDS Imaging System, 302nm*

GDST-1302 *GDS Touch Imaging System 302 nm*

**\$7,392**  
(List \$10,412)



WLV-12.19

You can now find Labnet on all major social media channels including Facebook, Twitter and YouTube



**FOLLOW US & LIKE US!**

**Labnet**  
Labnet International

33 Wood Ave. South, Iselin, NJ 08830  
(p) 732.417.0700

labnetinfo@corning.com  
www.labnetinternational.com

**Special offers valid January 1 - March 31, 2019**