

Prices and offers valid 1/1/21 through 3/31/21

RING IN THESE SAVINGS FOR THE NEW YEAR



LIQUID HANDLING ESSENTIALS

CONICAL TUBES



LPS Brand 15mL and 50mL Centrifuge Tubes



- Sterilized by gamma radiation
- Lip Seal Cap insures leak proof design
- Printed graduations and large writing area
- Polypropylene (PP) offer excellent chemical resistance and withstand temperature extremes of -40°C to 121°C. Maximum spin of 6,000 xg.
- Polystyrene (PS) offer glass-like clarity and tolerates aqueous solutions of mild bases or weak caids. Maximum spin of 1,000 xg.

Cat. #	Description	List	SALE
L226001	Centrifuge Tubes, 15mL, PP, STRL, grad, 50/rk, 500/CS	\$138.50	\$103.88
L226010	Centrifuge Tubes, 15mL, PP, STRL, grad, 50/bg, 500/CS	\$122.90	\$ 92.18
L226020	Centrifuge Tubes, 50mL, PP, STRL, grad, 25/rk, 500/CS	\$162.00	\$121.50
L226030	Centrifuge Tubes, 50mL, PP, STRL, grad, 25/bg, 500/CS	\$143.50	\$107.63

PCR® TUBES



PureAmp 0.2mL PCR® Tubes

- Certified RNase, DNase and PCR inhibitor FREE
- Ultra thin walls with uniform thickness
- Available with flat, domed or optical caps
- White tube option for better qPCR results



Cat. #	Description	List	SALE
L230105	PuerAmp qPCR 8-strip, 0.2mL w/opt clear strip caps, Wht 120	\$125.00	\$93.75
L230106	PuerAmp qPCR 8-strip, 0.1mL w/sep opt clear strip caps, 120	\$119.00	\$89.25
L230110	PureAmp indiv 0.2mL tubes w/dome caps, 1000/pk	\$40.00	\$30.00
L230109	PureAmp indiv 0.2mL tubes w/flat caps, 1000/pk	\$40.00	\$30.00

PCR® PLATES



UnSkirted 96-Well PCR® Plates

- Available in standard 200µL and low profile 100µL vol.
- Have high contrast alpha numeric printing
- Are thin wall and certified DNase, RNase FREE
- Have slightly raised rims for optimal sealing

Cat. #	Description	List	SALE
L226190	96-well PCR plate, 200µL, non-skirted, 100/UN	\$378.92	\$284.19

96-Well Semi-Skirted PCR® Plate

- Fit ABI® Cyclers and others
- Alphanumeric grid for easy sample identification
- RNase, DNase FREE

Cat. #	Description	List	SALE
L235900	96-well PCR plate, semi-skirted, fits ABI, 100/CS	\$398.42	\$298.82

384-Well Thin Wall PCR® Plates

- Designed for high throughput laboratory work
- Skirted to allow bar-coding on side for identification
- Each well has 40µl capacity
- Opaque plates ensure low levels of background fluorescents

Cat. #	Description	List	SALE
L234350	384-well, full skirt, A24, one notch, natural, 100/CS	\$523.98	\$332.98

PCR® FILM


PureAmp, PCR® Plate Sealing Film

- Compatible with all standard PCR Plates
- Available for PCR® and qPCR® Applications
- Certified RNase, DNase, and PCR inhibitor free
- Clear pressure sensitive adhesive insures evaporation rate of < 5%
- Packed 10 membranes per bag, 100 per pack

Cat. #	Description	List	SALE
T310011	PureAmp, sealing film, transparent, 10/pk, 100/cs	\$62.50	\$ 46.88
T310012	PureAmp, qPCR sealing film, pressure sensitive, 100/pk	\$149.00	\$111.75
T310014	PureAmp, Aluminum film, 10/pk, 100/CS	\$61.50	\$ 46.13



MICROCENTRIFUGE TUBES



Nest® Premium High Speed Locking Microtubes


- Available in 600µL, 1.5mL and 2.2mL with graduated markings
- Unique Safety Lock lid for secure sealing
- Rated to 30,000 xg, and certified RNase, DNase, DNA free

Cat. #	Description	List	SALE
605001	600µL Microtube with safe lock lid, graduated, frosted, 1000/PK	\$21.05	\$15.79
615001	1.5mL Microtube with safe lock lid, graduated, frosted, 500/PK	\$16.59	\$12.44
620011	2.2mL Microtube with safe lock lid, graduated, frosted, 500/PK	\$16.91	\$12.68

Micro Tube Storage Boxes

- Solid Polypropylene base with Clear Polycarbonate Cover
- 10x10 configuration for 600uL tubes and 9x9 for 1.5 / 2.2mL
- Temperature Range of -80 to +121°C, Freezer Safe
- Economical to use, 5 per pack

Cat. #	Description	List	SALE
613111	Microtube Storage Box for 1.5 & 2.2mL tubes, 9 x 9, 5 per pack	\$23.93	\$17.94



5mL MACRO SNAP CAP TUBES



5mL MacroTube™ Centrifuge Tubes

- Fit in standard 15mL centrifuge rotors or C2595 centrifuge
- Easy-open flip top with flat cap; provides tight leak-proof seal
- Syringe permeable cap
- 25,000xg; fully autoclavable
- RNase/DNase Free

Cat. #	Description	List	SALE
L225001	5.0mL MacroTube, EZ-Open, GRAD, NAT, NSTRL, 200/pk	\$39.50	\$29.63

5mL MacroTube™ Rack

- Holds 16 Macro tubes or 15mL Conical tubes in two tiers
- Non-skid rubber feet for stability
- Tubes are held securely in custom shaped wells
- Packed 2 per case

Cat. #	Description	List	SALE
T225920	MacroTube Rack for 5mL & 15mL Tubes, holds 16 2 per pack	\$49.50	\$37.13

CryoBox for the 5mL MacroTube™

- Polycarbonate freezer storage box will hold 25 tubes (5x5)
- Packed 5 boxes per pack
- Each box measures 5.25" x 5.25" x 2.25" OD

Cat. #	Description	List	SALE
T225810	CryoBox, Polycarbonate, 25 pl, 5mL tubes, 5 per pack	\$69.00	\$51.75

STERILE CULTURE TUBES

Sterile Culture Tubes with Two Position Caps



- 12 x 75mm and 17 x 100mm tubes with attached two position caps
- Packed individually, bags of 25 or bags of 125 per case
- Available with printed graduations and white marking area

Cat. #	Description	List	SALE
L288301	Culture Tube, Sterile, with cap 12 x75mm, PP, 25/bg,1000/cs	\$119.60	\$ 89.70
L288235	Culture Tube, Sterile, with cap 17x100mm, PP, 25/bg,1000/cs	\$209.00	\$156.75

DuoClick™ Screw-Cap Culture Tubes



- DuoClick™ tubes feature a unique, patent pending screw closure system.
- The advanced cap design includes two distinct positions for aerobic and anaerobic culturing
- Rotate the cap 1/4 turn to the first click to secure the cap in the aerobic position
- Rotate it another 1/2 turn past the second click for a leakproof, airtight seal for anaerobic culturing
- Available in 12 x 75 (5ml) or 17 x 100mm (14ml) sizes
- Sterile packed (20 tubes/bag or 50/foam rack), 500 per case

Cat. #	Description	List	SALE
L288330	DuoClick™ Screw-Cap Culture Tubes, 2-pos screw-cap, Strl, 17x100 bagged, 25/bag, 500/CS	\$194.00	\$145.50
L287300	DuoClick™ Screw-Cap Culture Tubes, 2-pos screw-cap, Strl, 12x75 bagged, 25/bag, 500/CS	\$135.00	\$101.25

RESERVOIRS

MTC Bio® ASPIR-8™ 5 & 10mL Reservoirs



- Made of white rigid polystyrene for added strength
- Molded in graduations with a tapered "V" bottom for easy pipetting
- Non-stick stacking design with pour off corners for easy use
- 5 & 10mL volumes are ideal for use with all 8-channel pipettors
- Available in non-sterile bulk pack, or sterile pack (individ. or 5 per)

Cat. #	Description	List	SALE
M170505	ASPIR-8 Reservoir, 5mL, Sterile, 5 per pk, 200/CS	\$104.00	\$78.00
M180105	Reagent Reservoir, 10 mL rigid styrene, Strl, 5/pk, 200/CS	\$104.00	\$78.00

Aquafil™ 25, 50 & 100mL Reservoirs



- Sharp -V- bottom eliminates waste
- Sterile pack, RNase, DNase, DNA Free and non-pyrogenic
- Unique stacking "lugs" for easy removal from stack
- Graduated with four pour off corners
- Available in 25mL, 50mL and 100mL capacity

Cat. #	Description	List	SALE
M255200	Reagent Reservoir, 25mL Sterile, 5/pack, 200/cs	\$122.68	\$ 92.01
M505200	Reagent Reservoir, 50mL Sterile, 5/pack, 200/cs	\$138.02	\$103.52
M100520	Reagent Reservoir, 100mL Sterile, 5/pack, 200/cs	\$143.12	\$107.34



Laboratory Products Sales Inc.

A NY State Certified Woman-Owned Small Business

1665 Buffalo Road • Rochester, N.Y. 14624 • 585-247-4720

TOLL FREE: 1-800-388-0166 • FAX: 585-247-6686 • www.LPSinc.com