

Prices and offers valid 10/1/20 through 12/31/20

# BACK TO BASICS

## Equipment and Supplies for Every Lab

### ACCUCOLD™ REFRIGERATORS AND FREEZERS

#### ACCUCOLD™ UNDER COUNTER FREEZER

- -25°C storage with extra thick insulation
- Slide-out drawers that help protect stored contents
- 3.5 cu.ft. 33.5"H x 23.5"W x 23.5"D
- 1 year full warranty, 5 yr. Compressor Warranty

Cat. #	Description	LIST	SALE PRICE
VT65ML	Summit 5 cu.ft. Freezer White, Under Counter, 115V	\$925.00	<b>\$695.00</b>



#### 5.5 CUBIC FT. REFRIGERATOR

- All-refrigerator, fully automatic defrost
- Front mounted lock
- Adjustable thermostat, no internal fans
- 33½"H x 23⅝"W x 23½"D

Cat. #	Description	LIST	SALE PRICE
FF-7L	Refrigerator w/lock, auto defrost, 5.5 cu.ft., white	\$970.00	<b>\$654.00</b>





### FREEZER RACKS AND STORAGE

#### UPRIGHT FREEZER RACKS FOR 2" & 3" BOXES

- Constructed of corrosion-resistant stainless steel
- Custom racks available on request • Locking rods available

Cat. #	Description	LIST	SALE PRICE
T396442	Upright rack for 16 x 2" boxes, 9¼"H x 22¼"D	\$102.20	<b>\$79.00</b>



#### FREEZER RACKS FOR 15 AND 50mL TUBE BOXES

- Racks Hold all 5¼" x 5¼" x 4¼" boxes containing tubes
- Constructed of corrosion resistant stainless steel
- Locking rods available

Cat. #	Description	LIST	SALE PRICE
T332320	Upright freezer rack for 15/50mL tube boxes, 3Wx2H	\$104.47	<b>\$83.58</b>

## COOL CELL® LX FREEZING CONTAINERS

- Provides alcohol-free freezing at the rate of -1°C/minute
- Holds 12 standard 1.0mL to 2.0mL cryogenic vials
- Radially symmetric for uniform vial freezing
- Available individually or as an assorted color 4 pack



Cat. #	Description	LIST	SALE PRICE
432002	CoolCell LX, Cell Freezing Container, Green 12 x 2mL, Ea	\$249.00	<b>\$186.75</b>



## 80 PLACE UNI-RACK® FOR PCR AND MICRO TUBES

- Polypropylene rack can be flipped to hold either PCR or 1.5/2mL tubes
- Racks can be anchored together or stacked for shelf storage
- Pack of 10 includes 2 ea. (Orange, Blue, Yellow, Green, Pink)

Cat. #	Description	LIST	SALE PRICE
T550080	Uni-Rack, 80 well flipper for PCR 0.5/1.5mL, ass't colors,10/pk	\$34.96	<b>\$26.22</b>

## ULT FREEZER BOX FOR 15 & 50mL TUBES

- Polypropylene box can take temperatures from -90 to 121°C
- Measures 144 x 144 x 127 mm, 102 mm without lid.
- Numerical index on grid
- 12 boxes per case



Cat. #	Description	LIST	SALE PRICE
T255701	Freezer rack, 36x15mL, w/lid, PP, green,12/pk	\$227.63	<b>\$159.00</b>



## TRUE NORTH® COOL CONTAINERS™

- Made from rigid molded polyurethane for durability
- Ideal for LN2, ice baths, ice-salt slurries and more
- Snug-fit lid keeps evaporation to a minimum
- Available in 4L bucket & 1,4,9 Liter color-coded pans

Cat. #	Description	LIST	SALE PRICE
HS28724B	True North™ Cool Container 4L Pan, Arctic Blue w/Lid, Ea	\$107.98	<b>\$72.00</b>

## PIPETTORS AND PIPET TIPS

### ALPHAPETTE PIPETTORS

- Lightweight, ergonomical and extremely durable
- Continuously adjustable to deliver the entire volume range
- Fully autoclavable, ideal for working under the hood
- Light and easy plunger action avoids hand fatigue
- Volume range from 2µL through 5mL
- Available in an economical starter pack with 4 pipets and a stand

Cat. #	Description	LIST	SALE PRICE
A-PCK	Alphapette Savings Bundle, incl 10,20,200,1mL pipets & stand	\$967.00	<b>\$595.00</b>
A2-PCK	Alphapette Savings Bundle, incl A2, 20, 200µl, 1mL pipet, stand	\$948.00	<b>\$595.00</b>



### DIAMOND RV-PETTE PRO REPEATER PIPET

- One-hand adjustment of tip volume for ease of use
- Accommodates 12 syringe tip sizes & dispenses up to 50mL
- Uses Universal RV-Pette Pro dispenser tips in sterile or non-sterile

Cat. #	Description	LIST	SALE PRICE
P39170	Diamond RV-Pette Pro, Repeat Pipettor, 120 volumes, EA	\$410.00	<b>\$295.00</b>



## LABNET BIOPETTE PLUS MULTI CHANNEL PIPETTERS

- Extremely accurate and fully autoclavable
- Reduced force tip ejector mechanism eases thumb stress
- Unique spring loaded cones makes positive tip seal • Available in 8 or 12 channel models



Cat. #	Description	LIST	SALE PRICE
P4508-200	BioPette Plus 20-200µL 8-channel pipettor	\$571.95	<b>\$423.00</b>
P4508-10	BioPette Plus 1-10µL 8-channel pipettor	\$571.95	<b>\$423.00</b>

## AVANT™ TIP STATION

- Ideal for multi-channel pipettor use • Double layer stack is ideal for use in hoods
- Stable base will not accidentally tip over • Precision fine point 10µL reach and 200µL pipet tips

Cat. #	Description	LIST	SALE PRICE
L122790	Avant Tip Station, sterile, 10µL reach, double stack, 960/un	\$34.10	<b>\$25.50</b>



## 10 L ULTRA MICRO-REACH TIPS

- Fit Gilson, Alhaphette, Finnpette, Eppendorf Research
- Low retention polymer ensures accurate pipetting
- Micro reach length fits to bottom of PCR® tubes
- Sterile Packed, 96 per box, 960 per pack




Cat. #	Description	LIST	SALE PRICE
L137642	Pipet Tip, 10µL reach, sterile, low ret. 96/rk, 960/PK	\$43.93	<b>\$34.35</b>

## TISSUE AND CELL CULTURE

### CELL CULTURE PLATES

- High clarity polystyrene; Thin wall design • Excellent bottom uniformity
- Non-reversible lid ensures low evaporation • Vacuum plasma TC treatment
- Easy peelable, individually wrapped package • Sterile; Non-pyrogenic

Cat. #	Description	LIST	SALE PRICE
703001	6-well Cell Culture Plate, Flat, TC, STRL, 1/pk, 50/CS	\$ 72.71	<b>\$54.54</b>
701001	96-well Cell Culture Plate, Flat, TC, STRL, 1/pk, 100/CS	\$147.46	<b>\$110.60</b>



### CELL CULTURE DISHES


- High clarity polystyrene • Flat transparent surface for distortion-free observation
- Vacuum plasma TC treatment; excellent cell adherence • PA/PE packaging; Stackable
- Vents for gas exchange • Sterile; Non-pyrogenic

Cat. #	Description	LIST	SALE PRICE
704004	Cell Culture Dish, 100x20mm, TC, STRL, 5/bag, 300/CS	\$109.88	<b>\$82.41</b>
705001	Cell Culture Dish, 60mm, TC, STRL, 20/pk, 500/CS	\$116.16	<b>\$87.12</b>

## NEST® CELL CULTURE FLASKS

- High clarity polystyrene • Vacuum plasma TC treatment for better cell attachment
- Option of plug seal cap or vent cap with 0.22µm hydrophobic filter
- Sterile; Non-pyrogenic

Cat. #	Description	LIST	SALE PRICE
708003	Cell Culture Flask, 75cm <sup>2</sup> , Vent, TC, STRL, 5/pk, 100/CS	\$150.70	<b>\$113.03</b>



## ERLENMEYER FLASKS, STERILE WITH SCREW CAPS

- PETG Flasks are individually wrapped with either plug or vent xaps
- Flasks are ideal for use in shaking incubator applications
- Vent caps have a 0.22µm filter for gas exchange without risk of contamination
- Available in 125, 250, 500 and 1000mL volumes.

Cat. #	Description	LIST	SALE PRICE
786011	3 Liter PC Erlenmeyer Flasks, Vent Cap, Sterile, 1/pk, 4/cs	\$244.43	<b>\$183.33</b>
787001	5 Liter PC Erlenmeyer Flasks, Seal Cap, Sterile, 1/pk, 4/cs	\$298.08	<b>\$223.56</b>

## PCR AND qPCR

### ULTRAFLUX™ IPCR STRIP TUBE WITH ATTACHED CAP

- Angularly attached caps prevent cap lids and hinges from interfering with each other
- Robust interwell linkages make the tubes rigid, decreasing the chance of spillage
- Certified free of detectable RNase DNase, DNA, PCR inhibitors, and tested pyrogen-free
- Packed 120 strips of 8 (960 tubes) per Pack



Cat. #	Description	LIST	SALE PRICE
L232470	UltraFlux 8-Strip PCR tube with attached flat cap, rigid strip	\$84.90	<b>\$63.68</b>



### NEST® 200 L PCR TUBE CAPS WITH FLAT CAPS

- Optically clear flat cap can be used for qPCR
- Thin walled tubes are Certified RNase, DNase, DNA Free
- Easy open leak proof caps reduce evaporation (<0.2%) • Packed 1000 tubes per pack


Cat. #	Description	LIST	SALE PRICE
401001	0.2mL PCR Tubes, Flat Caps, Clear, 1000/pk	\$41.86	<b>\$31.40</b>

## 96-WELL SEMI-SKIRTED PCR PLATE

- Fit ABI® Cyclers and others • Alphanumeric grid for easy sample identification
- RNase, DNase FREE

Cat. #	Description	LIST	SALE PRICE
L235900	96-well PCR Plate, semi-skirted, fits ABI, 100/CS	\$398.42	<b>\$298.82</b>





## THERMASEAL 2 PCR® SEALING FILM

- Heat resistant, recommended from -40°C to +120°C • Better adhesion; Reduced evaporation
- Longer end tabs for easier handling • certified DNase-, RNase-, and nucleic-acid-free

Cat. #	Description	LIST	SALE PRICE
T329020	ThermalSeal2 PCR film, 2mil thick, N-STRL, 100/PK	\$81.98	<b>\$61.49</b>

## PCRBio Dntp MIXES

- Each premixed aqueous solution contains dATP, dCTP, dGTP, dTTP
- Available in 10mM or 25mM volumes
- Ultra Pure (>99%) and remains stable after repeated freeze/thaw cycles
- Ideal for all PCR and qPCR applications and sequencing

Cat. #	Description	LIST	SALE PRICE
PB10.72-10	dNTP Mix 25mM each (100mM total) 10µmol (1 ml)	\$130.00	<b>\$104.00</b>



## qPCRBio SYGREEN MIXES

- Non-PCR inhibiting intercalating dye • Rapid extension rate for early Ct values
- Market leading sensitivity • Available with clear or easy to see Blue mix
- Choose Hi-ROX, Lo-ROX, Separate-ROX or with Fluorescein

Cat. #	Description	LIST	SALE PRICE
PB20.15-20	SyGreen Blue Mix Lo-ROX 2000 x 20µL rxns	\$810.00	<b>\$607.50</b>

## DISINFECTANT CLEANERS AND WIPES

### EASYTASK® REFILLABLE WIPE SYSTEM

- Solvent resistant durable wipes in a perforated centerfeed roll
- Refillable GrabBox™ holds 2 quarts (2.2 liters) of disinfectant
- Resealable bag keeps wipes moist • Wipe refill rolls come 6 per CS

Cat. #	Description	LIST	SALE PRICE
6100255	EasyTask GrabBox Wipes, Centerfeed Roll, 275/roll, Box	\$ 50.00	<b>\$38.00</b>
6100257	EasyTask GrabBox Refill Wipes, 275/roll, 6 rolls/CS	\$230.00	<b>\$172.50</b>



### ECOLAB® PEROXIDE MULTI SURFACE CLEANER & DISINFECTANT

- Broad surface disinfectant, kills bacteria and virus (1:32 dilution)
- Ideal for use in schools and colleges, clinics and animal life science labs
- Available in a two gallon concentrate, No PPE required when using diluted

Cat. #	Description	LIST	SALE PRICE
6100693	Ecolab Peroxide Surface Cleaner & Disinfectant, 2 Gal	\$133.20	<b>\$99.90</b>



## Laboratory Products Sales Inc.

*A NY State Certified Woman-Owned Small Business*

1665 Buffalo Road • Rochester, N.Y. 14624 • 585-247-4720

TOLL FREE: 1-800-388-0166 • FAX: 585-247-6686 • www.LPSinc.com