

Prices and offers valid from 4/1/2020 through 6/30/2020



**Labnet**  
Labnet International, Inc.

# SPRING IS HERE! LAB ESSENTIALS SAVE UP TO 75%

## MultiGene Gradient Thermal Cycler

- Produces consistent, reliable PCR amplifications while making gradient technology affordable
- The MultiGene Gradient Thermal Cycler gradient is spread horizontally allowing twelve temperatures within a 24°C gradient
- 96 x 0.2mL block and an adjustable heated lid to fit a variety of PCR plates and PCR tubes
- Also available with a 60 x 0.5mL tube block model

**SAVE \$665**



Cat. #	Description	List	SALE
TC9610	MultiGene Gradient Thermal Cycler	\$3,330.74	<b>\$2,665.00</b>

## Enduro Semi-Dry Blotter

- For Western, Southern and Northern blots
- Rapid transfer time; Blotting can be done in 15-30 minutes
- Uniform head dispersion; Variable gel thickness
- Dimensions: 12.8x9.8x2.2" (32.5x25x5.5cm)
- 20mL Buffer volume
- 4 Blots: 8x8.5cm - 2 Blots: 16x8.5cm - 1 Blot: 16x17.5cm
- For total blotting solutions, use with ENDURO™ 250V, 3000mA high current power supply (E-0203/not included)



Cat. #	Description	List	SALE
E2020-SDB	Enduro Semi-Dry Blotter	\$1,228.84	<b>\$978.00</b>

## Enduro Power Supplies

- Dependable Performance, Intelligent control
- Constant Voltage, Current or Wattage
- Three high performance models for the most demanding jobs

Cat. #	Description	List	SALE
E-0303	Enduro 300V Power Supply, 2-300V, 4 to 500mA, 4 outlets, 120V	\$498.66	<b>\$403.00</b>
E-0203	Enduro 250V Power Supply, 10-3000V, mA programmable, 4 outputs, 250V	\$962.37	<b>\$702.24</b>

## Labnet Mini

### Microcentrifuge

- Includes 8-place rotor for 1.5 & 2.0mL tubes
- Spins at 6000 rpm and 2000xg
- Available with Blue, Red, Green or Purple accented lids
- Electronic braking with dual lid switch



Cat. #	Description	List	SALE
C-1601-B	Labnet Mini Microcentrifuge, Blue lid	\$195.00	<b>\$156.00</b>
C-1601-P	Labnet Mini Microcentrifuge, Purple lid	\$195.00	<b>\$156.00</b>
C-1601-R	Labnet Mini Microcentrifuge, Red lid	\$195.00	<b>\$156.00</b>
C-1601-G	Labnet Mini Microcentrifuge, Gray lid	\$195.00	<b>\$156.00</b>

## Labnet BioPette™ Plus Single Channel Pipettors

- Superior accuracy and reproducibility
- Adjustable volume with thumbwheel or push-button
- Fully autoclavable with integral volume lock



Cat. #	Description	List	SALE
P3942-2	BioPette™ Plus Pipettor, 0.1-2.0µL	\$240.00	<b>\$180.00</b>
P3942-10	BioPette™ Plus Pipettor, 0.5-10µL	\$240.00	<b>\$180.00</b>
P3942-20	BioPette™ Plus Pipettor, 2-20µL	\$226.00	<b>\$169.50</b>
P3942-100	BioPette™ Plus Pipettor, 10-100µL	\$226.00	<b>\$169.50</b>
P3942-200	BioPette™ Plus Pipettor, 20-200µL	\$226.00	<b>\$169.50</b>
P3942-1000	BioPette™ Plus Pipettor, 100-1000µL	\$226.00	<b>\$169.50</b>
P3942-5000	BioPette™ Plus Pipettor, 1000-5000µL	\$228.26	<b>\$171.20</b>
P3942-10000	BioPette™ Plus Pipettor, 1000-10000µL	\$264.01	<b>\$198.01</b>

## Labnet BioPette™ Plus Starter Kit

- 4 Pipettes (0.5-10 µL, 20-200 µL and 100-1000µL)
- 3 Racks of BioFree Tips
- 6-place Carousel Stand

Cat. #	Description	List	SALE
P3942-SK4	BioPette™ Plus Starter Kit	\$818.40	<b>\$613.80</b>



FOLLOW US ON:



**Laboratory Products Sales Inc.**

*A NY State Certified Woman-Owned Small Business*

1665 Buffalo Road • Rochester, N.Y. 14624 • 585-247-4720

TOLL FREE: 1-800-388-0166 • FAX: 585-247-6686 • www.LPSinc.com