

### Z287-A Microcentrifuge PLUS

- Bundle includes 2 rotors  
18 x 1.5 / 2.0 ml and 2 x microplate
- Up to **14,000 rpm** and **16,058 xg**
- Indicator Light, easily view cycle status

#### Bundle Includes



Z287-02MP  
2 x microplate  
3,800 rpm / 1,453 xg



Z287-1820  
18 x 1.5 / 2.0 ml  
14,000 rpm / 15,994 xg

#### Also Available



6 x 50 ml  
6,000 rpm / 4,427 xg



24 x 1.5 / 2.0 ml  
12,000 rpm / 13,683 xg



Unbeatable value!  
**\$3795**

Compact, only 14in. wide

#### BUNDLE ORDERING INFORMATION

	Reg.	Bundle
Z287-BND-MB	Z287-A Micro bundle	\$5530 <b>\$3795</b>
(incl. Z287-1820 and Z287-02MP)		

### Benchtop Universal Centrifuge Bundles



As low as **\$6,750**

#### Z326 Series

- Max. swing-out capacity:  
8 x 50ml or 12 x 15ml
- Truly universal, 17 rotor options



As low as **\$8,200**

#### Z366 Series

- Max. swing-out capacity:  
16 x 50ml or 40 x 15ml
- High speed rotors up to **24,000 xg**



As low as **\$9,800**

#### Z446 Series

- Max. swing-out capacity:  
28 x 50ml or 68 x 15ml
- Microplate capacity:  
20 x microplates or 8 x PCR Plates

#### BUNDLE ORDERING INFORMATION

		Reg.	Bundle
Z326-BND	Z326 Ambient Swing Out Bundle	<del>\$7,618</del>	<b>\$6,750</b>
Z326-K-BND	Z326-K Refrigerated Swing Out Bundle	<del>\$10,136</del>	<b>\$9,100</b>
Z366-BND	Z366 Ambient Swing Out Bundle	<del>\$9,087</del>	<b>\$8,200</b>
Z366-K-BND	Z366-K Refrigerated Swing Out Bundle	<del>\$12,136</del>	<b>\$10,600</b>
Z446-BND	Z446 Ambient Swing Out Bundle	<del>\$10,861</del>	<b>\$9,900</b>
Z446-K-BND	Z446-K Refrigerated Swing Out Bundle	<del>\$14,121</del>	<b>\$13,000</b>
Z446-MPBND	Z446 Ambient Microplate Bundle	<del>\$10,903</del>	<b>\$9,800</b>
Z446-K-MPBND	Z446-K Refrigerated Microplate Bundle	<del>\$14,183</del>	<b>\$12,500</b>



\*Additional bundles available upon request.  
To order 230V voltage, add "E" to part number.

### Z306 Universal Centrifuge



**\$5,300**

- Includes swing-out rotor for up to **8 x 50ml** or **12 x 15ml**
- Universal, 9 rotor options
- Compact, less than 14 inches wide



Swing-out bundle

#### BUNDLE ORDERING INFORMATION

Z306-BND	Z306 Universal centrifuge with swing-out rotor	Reg. <del>\$5,877</del>	<b>Bundle \$5,300</b>
----------	--	-------------------------	-----------------------

Reg. **Bundle**

### Z216-MK High Speed Refrigerated Microcentrifuge

**REFRIGERATED**



**As low as \$5,600**

- **2** bundle options!
- Compact, only 11 inches wide
- High speed, up to **21,631xg**
- **44** place rotor option
- **COMBI-rotor™** option with Quick-Seal™ lid

**NEW** Hermle COMBI-rotor™ with Quick-Seal lid



Holds 24 x 1.5/2.0ml and 4-PCR strips

**Quick-Seal™ lid**  
Lid opens and locks with only a quarter turn



Now included in all Z216 series bundles!

#### BUNDLE ORDERING INFORMATION

Z216-K-BND44	Z216-MK Refrigerated Microcentrifuge w/ 44 place rotor	Reg. <del>\$6,569</del>	<b>Bundle \$5,600</b>
Z216-K-CMB	Z216-MK Refrigerated Microcentrifuge w/ COMBI-Rotor™	Reg. <del>\$7,120</del>	<b>Bundle \$5,800</b>

Reg. **Bundle**

\*Ambient versions also available contact your distributor for pricing

### Z36 HK Super Speed Centrifuge

\*Additional bundles available upon request. To order 230V voltage, add "-E" to part number.



**As low as \$14,000**

- Speeds up to **65,395 xg!**
- Max. swing-out capacity: **12 x 50ml** or **28 x 15ml**
- Refrigerated, **-20°C to 40°C**



6 x 250ml rotor  
10,000 rpm / 15,650 xg



6 x 50ml rotor  
21,000 rpm / 41,410 xg



12 x 1.5/2.0 ml rotor  
30,000 rpm / 65,395 xg

Purchase a Z36-HK BUNDLE and receive a **FREE BenchPACK!™**

#### BUNDLE ORDERING INFORMATION

Z36-BND250	Z36-HK Centrifuge with 6 x 250ml rotor	Reg. <del>\$18,652</del>	<b>Bundle \$15,000</b>
Z36-BND50	Z36-HK Centrifuge with 6 x 50ml rotor	Reg. <del>\$16,923</del>	<b>Bundle \$14,000</b>
Z36-BND20	Z36-HK Centrifuge with 12 x 1.5/2.0ml rotor	Reg. <del>\$16,131</del>	<b>Bundle \$14,000</b>

Reg. **Bundle**

#### INCLUDES:

- 1 **myFuge™ 12** Mini Centrifuge
- 1 **BenchMixer™ V2** Vortex Mixer
- 1 **myTemp™ Mini** Incubator
- 1 **Accuris™ Mini** Balance
- 1 **Backpack**

**Value \$1,500**

