

## Z206-A Compact Centrifuge



- Includes **6 x 50ml** rotor
- Compact, only 11 inches wide
- Optional rotors, for tubes from **0.2ml to 50ml**
- Bundle includes **FREE** ProPette™ Electronic pipette controller with QuickStand™

**\$2255**

6,000 rpm / 4,427 xg

Max. Speed: 6,000 rpm  
Max. g-Force: 4,427 xg  
Max. Capacity: 6 x 50ml  
12 x 15ml  
Dimensions: 11 x 14.5 x 10"

**FREE!**

**QUICKSTAND™**

\*with purchase of Z206-A centrifuge bundle

### BUNDLE ORDERING INFORMATION

Z206-BND	Z206-A Compact centrifuge with 6 x 50ml rotor	<del>\$2884</del>	<b>\$2255</b>
P6080	ProPette™ Pipette Controller, w/ QuickStand™	<del>\$279</del>	<b>FREE</b>
Z206-50-A15	Optional adapters for 15ml tubes, pack of 6	\$118	

## Z306 Universal Centrifuge



**Microplate bundle**

Includes **2 x 3 Microplate Rotor**  
4,500 rpm / 2,716 xg

**\$5150**

- Includes 2 x 3 microplate rotor
- Universal, 9 rotor options
- Compact, less than 14 inches wide
- Optional, swing-out rotor for **8 x 50ml** or **12 x 15ml**

Max. Speed: 14,000 rpm  
Max. g-Force: 18,624 xg  
Max. Capacity: 8 x 50ml  
12 x 15ml  
Dimensions: 14 x 18.6 x 13"

Universal, chamber slide and 35/60mm dish carrier for microplate rotors



Item: Z326-MPCA \$189, pack of 2

**Swing-out bundle**



Optional swing-out rotor (Included with Z306-BND)

### BUNDLE ORDERING INFORMATION

Z306-BND-MP	Z306 Universal centrifuge with microplate rotor	<del>\$5734</del>	<b>\$5150</b>
Z306-BND	Z306 Universal centrifuge with swing-out rotor	<del>\$5823</del>	<b>\$5050</b>

\*Limited Time offer, valid through June 30th 2019

## Z216-M High Speed Microcentrifuges

- Includes **44 place 1.5/2.0ml rotor**
- Available with and without **refrigeration**
- Compact, only 11 inches wide
- High speed (21,379xg with 44 place rotor)
- Optional rotors for 5ml tubes down to 0.2ml tubes/strips



Z216-M

**\$3050**



Z216-MK  
Refrigerated

**\$5395**

15,000 rpm / 21,379 xg

Max. Speed: 15,000rpm  
Max. g-Force: 21,379 xg  
Max. Capacity: 44 x 1.5/2.0ml  
Temp. Range: -20 to 40°C (MK)  
Dimensions: 11 x 15 x 11" (M)  
11 x 21 x 11" (MK)

### BUNDLE ORDERING INFORMATION

Z216-BND-44	Z216-M Microcentrifuge with 44 place rotor	<del>\$3609</del>	<b>\$3050</b>
Z216-K-BND-44	Z216-MK Refrigerated Microcentrifuge with 44 place rotor	<del>\$6409</del>	<b>\$5395</b>
R1000	CoolCube™ Microtube and PCR Plate Cooler	<del>\$114</del>	<b>FREE</b>

**Special OFFER**

Item: R1000



\*with purchase of Z216 centrifuge bundle

**FREE!**

\*Limited Time offer, valid through June 30th 2019

### 4 x 145ml Capacity **Z326 Series** As low as **\$6725!**

- Swing-out capacity: **8 x 50ml** or **12 x 15ml**
- Truly universal, 17 rotor options

Max. Speed: 18,000 rpm      Max. g-Force: 17,317xg  
 Max. Capacity: 4 x 145ml      Temp. Range: -20 to 40°C (-K)  
 Dimensions: 16 x 19 x 14"      16 x 27.5 x 14" (Z326-K)

#### BUNDLE ORDERING INFORMATION

Z326-BND      Z326 Universal centrifuge BUNDLE\*      ~~\$7434~~      **\$6725**  
 Z326-K-BND      Z326-K Refrigerated Centrifuge BUNDLE\*      ~~\$9889~~      **\$9050**



\*Bundle Includes: Centrifuge, 4 x 100ml swing out rotor, 3x15ml bucket (pack of 2) and 2x50ml bucket (pack of 2)

### As low as **\$8095!** **Z366 Series** 4 x 250ml Capacity



- Swing-out capacity: **16 x 50ml** or **40 x 15ml**
- High speed rotors available up to 24,000 x g

Max. Speed: 16,000 rpm      Max. g-Force: 24,325xg  
 Max. Capacity: 4 x 250ml      Temp. Range: -20 to 40°C (-K)  
 Dimensions: 17 x 20 x 14"      28 x 20 x 14" (-K)

#### BUNDLE ORDERING INFORMATION

Z366-BND      Z366 Universal Centrifuge BUNDLE\*      \$8866      **\$8095**  
 Z366-K-BND      Z366-K Refrigerated Centrifuge BUNDLE\*      ~~\$11867~~      **\$10700**

\*Bundle Includes: Centrifuge, 4 x 250ml swing out rotor, 10x15ml bucket (pack of 2) and 4x50ml bucket (pack of 2)

### 4 x 750ml Capacity **Z446 Series** As low as **\$9685!**

- Swing-out capacity: **28 x 50ml** or **68 x 15ml**
- Microplate capacity: **12 x microplates** or **8 x PCR Plates**

Max. Speed: 16,000 rpm      Max. g-Force: 24,325xg  
 Max. Capacity: 4 x 750ml      Temp. Range: -20 to 40°C (-K)  
 Dimensions: 21 x 26.5 x 15.5"      29 x 26.5 x 15.5" (-K)

#### BUNDLE ORDERING INFORMATION

Z446-BND      Z446 Tissue Culture (TC) BUNDLE\*      \$10597      **\$9755**  
 Z446-K-BND      Z446-K Refrigerated Tissue Culture (TC) BUNDLE\*      \$13797      **\$12670**  
 Z446-MPBND      Z446 Microplate (MP) BUNDLE\*      \$10638      **\$9685**  
 Z446-K-MPBND      Z446-K Refrigerated Microplate (MP) BUNDLE\*      \$13838      **\$12650**

**TC Bundle Includes:** Centrifuge, swing-out rotor, 17x15ml bucket (pk of 2) & 7x50ml bucket (pk of 2)  
**MP Bundle Includes:** Centrifuge, 4 x 750ml swing-out rotor, 4 microplate carriers



Tissue culture bundle

\* Lids not included

Microplate bundle

Z446-750  
Swing out rotor      4,500 rpm / 4,346 xg

#### Also Available

**Z446-MP-BUK\***  
Rectangular bucket,  
for microplates (sealable)



\*Lid not included

**Z446-MP-STK**  
PCR Plate Stacking Adapter,  
for Z446-750-MP



Stacks PCR plates for unmatched capacity!