

LAB VALUES

from
Labnet International

FALL 2018
LIQUID HANDLING

Special offers valid September 1- December 31, 2018

ACCURACY AND PRECISION

Volume Range	Accuracy	Precision
0.1-2.0 µl	±12 to ±1.5%	≤6.0 to ≤0.7%
0.5-10 µl	±4.0 to ±0.5%	≤2.8 to ≤0.4%
2-20 µl	±3.0 to ±0.8%	≤1.5 to ≤0.3%
10-100 µl	±1.6 to ±0.8%	≤0.8 to ≤0.2%
20-200 µl	±1.2 to ±0.6%	≤0.6 to ≤0.2%
100-1,000 µl	±1.6 to ±0.6%	≤0.40 to ≤0.15%
1,000-5,000 µl	±0.6 to ±0.5%	≤0.25 to ≤0.15%
1,000-10,000 µl	±2.5 to ±0.5%	≤0.6 to ≤0.2%

SINGLE CHANNEL PIPETTE

Ordering Information

Cat No.	Volume Range	Tip	List	Special
P3942-2	0.1-2.0µl	10µl	\$240	\$192
P3942-10	0.5-10µl	10µl	\$240	\$192
P3942-20	2-20µl	200µl	\$226	\$181
P3942-100	10-100µl	200µl	\$226	\$181
P3942-200	20-200µl	200µl	\$226	\$181
P3942-1000	100-1,000µl	1,000µl	\$226	\$181
P3942-5000	1,000-5,000µl	5,000µl	\$226	\$181
P3942-10000	1,000-10,000µl	10,000µl	\$264	\$211

8 CHANNEL PIPETTE

Cat No.	Volume Range	Tip	List	Special
P4808-10	1-10µl	10µl	\$572	\$458
P4808-50	5-50µl	200µl	\$572	\$458
P4808-200	20-200µl	200µl	\$572	\$458
P4808-300	50-300µl	200/300µl	\$572	\$458

12 CHANNEL PIPETTE

Cat No.	Volume Range	Tip	List	Special
P4812-10	1-10µl	10µl	\$682	\$546
P4812-50	5-50µl	200µl	\$682	\$546
P4812-200	20-200µl	200µl	\$682	\$546
P4812-300	50-300µl	200/300µl	\$682	\$546



SAVE 20% or more

See Page 4
for
Clearance
Items

Labmax™ Bottletop Dispensers



D5370-10	2-10 mL	\$280
D5370-10HF*	2-10 mL for HF	\$788
D5370-2	0.3-2.5 mL	\$280
D5370-25	5-25 mL	\$377
D5370-25HF*	5-25 mL for HF	\$1,046
D5370-5	1-5 mL	\$280
D5370-50	10-50 mL	\$406

*HF=PTFE piston and cylinder version for HF and other aggressive solvents

FLV-2.18

BioFree™ Pipette Tips

.5-10 uL tips		
P1000-10	Pk 1,000	\$388
P1000-10-R	10 racks of 96	\$222
200 uL tips		
P1000-200	Pk 1,000	\$388
P1000-200-R	10 racks of 96	\$222
300 uL tips		
P1000-300-R	10 racks of 96	\$222
1,000 uL tips		
P1000-1000-R	10 racks of 96	\$250
1,200 uL tips		
P1000-1200-R-F	10 racks of 96 filtered	\$230
10,000 uL tips		
K943	Bulk pack of 100	\$36



FLV-3.18

Labpette R™ Repeating Pipette

- Repetitive pipetting from single filling of disposable tip
- Just select tip & lock in place
- Chart located on body showing volume dispensed for each tip size
- Best used with BioFree Combi-Syringe non sterile or sterile dispenser tips

\$274



Ordering Information
P3000 **LabPette R Repeating Pipette**
1 µl - 5,000 µl

FLV-4.18

BioFree™ Combi-Syringe Tips

Tip	Cat No. Nonsterile	Price	CAT No. Sterile	Price
0.05 mL	P3505	\$170	P3505-S	\$191
0.5 mL	P3512	\$81	P3512-S	\$100
1.25 mL	P3515	\$80	P3515-S	\$100
2.5 mL	P3518	\$80	P3518-S	\$100
5.0 mL	P3521	\$80		
12.5 mL	P3524	\$80	P3524-S	\$100
25.0 mL	P3525*	\$109		
50.0 ml	P3527*	\$113	P3527-S	\$130
	P3535	Combination pk, 20 ea of 0.5, 1.25, 5.0 and 12.5 mL		\$80
	P3527-A	Adapter for 25 & 50 mL		\$7

* Four adapters included


FLV-5.18

FastPette™ V2 & FastPette™ Pro Pipette Controllers

- Ergonomically located switches
- Fully charged can be used continuously for 1,000 pipetting cycles
- Designed for use with glass or plastic serological pipettes from 0.5 to 100mL range
- NiMH batteries when fully charged - 8 hours continuous use

\$239

\$300



FLV-6.18

Ordering Information

P2000 **Fastpette V2** Pipette Controller
P2002 **Fastpette Pro** Pipette Controller

BioPette™ Plus Starter Kit

FREE
5th single channel pipette
with purchase of starter kit
\$654

Ordering Information

P3942-SK4 **BioPette Plus** Four pack starter kit includes: 4 pipettes (0.5-10 µl, 2-20 µl and 100-1000 µl) 3 racks of BioFree tips and a 6 place carousel stand

FLV-7.18

Excel™ Electronic Pipettes

	<i>LIST</i>	<i>SPECIAL</i>
Single channel 0.5-1200 µl	\$346	\$277
8-channel 0.5-200 µl	\$706	\$565
8-channel 1200 µl	\$800	\$640
12-channel 20-200 µl	\$824	\$659

- Simple, precise volume setting
- Fully motorized piston drive
- Multifunctional with 6 operating modes
- Virtually eliminates pipetting force
- Rechargeable lithium-ion battery



Ordering Information

P3600L-10	0.5 - 10 µl single channel
P3600L-20	2 - 20 µl single channel
P3600L-200	10 - 200 µl single channel
P3600L-1200	100 - 1200 µl single channel
P3608L-10	0.5 - 10 µl eight channel
P3608L-20	2 - 20 µl eight channel
P3608L-200	10 - 200 µl eight channel
P3608L-1200	100 - 1200 µl eight channel
P3612L-20	2 - 20 µl twelve channel
P3612L-200	10 - 200 µl twelve channel

FLV-8.18

Labpette FX™ Fixed Volume Pipettes

- Fully autoclavable
- Ergonomically designed handle and finger rest
- Unique "twist-eject" tip removal system reduces stress

\$80



Ordering Information

<u>200µL tip</u>		<u>1,000µL tip</u>	
P3900-5	P3900-70	P3900-200	P3900-500
P3900-10	P3900-75	P3900-220	P3900-600
P3900-20	P3900-80	P3900-250	P3900-700
P3900-25	P3900-90	P3900-300	P3900-800
P3900-30	P3900-100	P3900-400	P3900-900
P3900-40	P3900-120	P3900-450	P3900-1000
P3900-50	P3900-150		



FLV-9.18

High Performance Benchtop Centrifuge

FLV-10.18

- Ambient and refrigerated models
- Capacity up to 400 mL with swing out & fixed angle rotor
- Speed range 90-18,000 rpm/24,270 x g



C0336R

\$2,861
Ambient

\$4,959
Refrigerated

C0336

Ordering Information

C0336 High Performance Bench Top Centrifuge, ambient
C0336R High Performance Bench Top Centrifuge, refrigerated

Mini Power Supply

FLV-12.18

- Timer 1 m-99 h 59 m or continuous
- 2.5 lbs
- Small footprint 6 x 8 x 3"
- Optional Power Stack available

\$322



Ordering Information

E0304 Mini Power Supply

Labnet Mini Centrifuge FLV-11.18



- Quick release rotor system
- Electronic braking
- Compact, rugged design
- Dual lid switch
- Rotors accommodate 8 microtubes or 4 PCR strip tubes

\$156

Ordering Information

C1601* Mini Centrifuge, with 8 place tube rotor and 4 place PCR strip tube rotor & adapters for 0.5 mL tubes, 0.2 mL tubes and 0.25 mL tubes
*add-B (blue), -G (gray), -P (purple), -R (red), for lid color

Enduro™ GDS/GDST* FLV-13.18

*Gel Documentation System

- Plug & play installation < 15 min
- Large 20 x 24 field of view
- 4 click image acquisition
- GDST-blue light compatible with most DNA safe dyes



\$4,995

\$6,470

Laptop not included

Ordering Information

GDS-1302 GDS Imaging System, 302nm
GDST-1302 GDS Touch Imaging System, 302 nm

LAST CHANCE

25% discount on discontinued units below

Order #	Speed Range	Maximum Capacity	Price	
C1301* *add -B (blue) -P (purple)	6,000 rpm	2 x 0.2 mL strips or 6 x 1.5/2.0 mL tubes	\$146	FLV-14.18
C0206-A	200-6,000 rpm	12 x 15 mL or 6 x 50 mL	\$1,478	FLV-15.18
C0306	10-13,500 rpm	4 x 100 mL	\$2,638	FLV-16.18
C0326	200-18,000 rpm	4 x 100 mL	\$3,878	FLV-17.18
C0446	200-16,000 rpm	4 x 750 mL	\$5,134	FLV-18.18

You can now find Labnet on all major social media channels including Facebook, Twitter and YouTube



FOLLOW US & LIKE US!

Labnet

Labnet International

33 Wood Ave. South, Iselin, NJ 08830
(p) 732.417.0700
labnetinfo@corning.com
www.labnetinternational.com