

2025 LabMart Research Essentials Summer into Fall Savings



QuickSILVER™ Powdered Buffer Packs

- Each box makes 50 or 100 liters of 1X buffer
- Quality and consistency, guaranteed
- 80% less shipping weight than concentrates

EB1000	QuickSilver TAE Buffer Powder, 100 pouches	\$277
EB1100	QuickSilver Tris-Glycine-SDS Buffer Powder, 50 pouches	\$197
EB1200	QuickSilver PBS Buffer Powder, 100 pouches	\$177
EB1201	QuickSilver PBS-T Buffer Powder, 100 pouches	\$177
EB1202	QuickSilver TBS Buffer Powder, 100 pouches	\$247
EB1203	QuickSilver TBS-T Buffer Powder, 100 pouches	\$247



**Stop wasting valuable time...
Prepare liter quantities of
fresh 1X buffer in minutes**

myVolt™ Touch Electrophoresis Power Supply



- Touch screen control
- Dual orientation, auto-rotating display
- High performance specifications:
 - Max. Voltage: 300V
 - Max. Amperage: 400mA
 - Max. Wattage: 100W

E2301 myVolt Touch Power Supply **\$679**

Accuris UV Transilluminator

- Large, 21x26cm viewing surface
- Powerful 14 watt bulbs
- Variable intensity
- 302nm midrange UV wavelength
- Optimal for EtBr, SmartGlow™ and SYBR stains
- Hinged, UV-blocking cover



E3200 UV Transilluminator, Variable Intensity **\$1,250**

Accuris Analytical Balances

- Readability of 0.0001g
- Available with Quick-Cal™ internal calibration feature
- Multiple mode options (standard weighing, count mode, & percent)



Analytical Balances			Compact Balances		
W3100-120	120g x 0.0001g	\$2,164	W3300-120	120g x 0.01g	\$345
W3100A-120	120g x 0.0001g w/ QuickCal™	\$2,296	W3300-500	500g x 0.01g	\$440
W3100-210	210g x 0.0001g	\$2,513	W3300-1200	1200g x 0.01g	\$447
W3100A-210*	210g x 0.0001g w/ QuickCal™	\$2,663			

Free!*
With purchase
of W3100A-210

Accuris Compact Balances

- Readability of 0.1g to 0.01g
- Easy in-lab calibration feature
- RS232 port for data output to a computer or printer

***FREE!**
**Compact
Balance**



Use promo code: **LM225** With purchase of W3100A-210

Solution Reservoirs Standard 12-Channel Width

The LabMart®
Serving Science Since 1850

Solution Reservoirs Reduced 8-Channel Width

- Graduated inner walls
- Micro-Trough U-shaped bottoms
- Pouring spouts on all 4 corners
- RNase, DNase, & pyrogen free
- Pack of 5 protective lids with standard 25 and 50mL reservoirs



25mL - Shown with lid 50mL - Shown with lid 100mL

Description	25mL* Standard	50mL* Standard	50mL Divided	100mL Standard	QTY	Price
Non-Sterile Bulk Pack	ER-025NS	ER-050NS	ER-051NS	ER-100NS	100	\$64.36
Sterile Individually Wrapped	ER-025S-1	ER-050S-1	ER-051S-1	ER-100S-1	100	\$134.91
Sterile Packs of 5	ER-025S-5	ER-050S-5	ER-051S-5	ER-100S-5	200	\$134.91
Lids - Sterile Packs of 5	ER-025LID	ER-050LID	N/A	N/A	200	\$93.42

*5 FREE Lids included with every purchase of Standard 25mL and Standard 50mL reservoirs.

75mm Inner Chamber



For 8-Channel Pipettes

- Specifically designed for 8-channel pipettes
- 75mm inner chamber
- Conserve valuable reagents

Description	5mL 8-Channel	10mL 8-Channel	25mL 8-Channel	QTY	Price
Non-Sterile Bulk Pack	ER-005NS	ER-010NS	ER-025NS	300*	\$122.98
Sterile Individually Wrapped	ER-005S-1	ER-010S-1	ER-025S-1	100	\$122.98
Sterile Packs of 5	ER-005S-5	ER-010S-5	ER-025S-5	200	\$122.98

*ER-025NS has quantity of 100/pack and price of \$61.49

dot scientific inc.

Cryogenic Vials Internal or External Thread

- Withstands temperatures as low as -196°C
- Self-Standing star foot design
- Leakproof silicone O-ring
- Sterile; Resealable bags of 50
- Paper/plastic packaging



Externally Threaded Cryogenic Vials



Internally Threaded Cryogenic Vials

External	Internal	Description	Qty	Price
V3801	V4801	1mL Cryovial, Sterile, Silicone O-Ring, Graduated, Standing	500	\$164.16
V3802	V4802	2mL Cryovial, Sterile, Silicone O-Ring, Graduated, Standing	500	\$172.94
V3803	V4803	3mL Cryovial, Sterile, Silicone O-Ring, Graduated, Standing	500	\$219.93
V4805	V4805	5mL Cryovial, Sterile, Silicone O-Ring, Graduated, Standing	500	\$248.13
Cryogenic Vial Accessories				
V5809-A*		Cap Inserts for Cryogenic Vials, Assorted Colors	500	\$20.89
V5920		Locking rack for cryogenic vials, 50-place (5x10)	4	\$56.86

Serological Pipettes


- Positive and negative graduations
- Sterile, individually wrapped
- Cotton plugged
- Standardized color coding
- Paper/plastic packaging



Cat No.	Description	Qty	Price
P8270-01	1mL Pipette, Yellow Code, Sterile	500	\$247.57
P8270-02	2mL Pipette, Green Code, Sterile	500	\$247.57
P8270-05	5mL Pipette, Blue Code, Sterile	200	\$247.57
P8270-10	10mL Pipette, Purple Code, Sterile	200	\$247.57
P8270-25	25mL Pipette, Red Code, Sterile	150	\$247.57
P8270-50	50mL Pipette, Black Code, Sterile	150	\$247.57

SureStrain™ Cell Strainers Mesh Strainers in a Polypropylene Frame

- 8 pore sizes for greater cell sorting precision
- Designed for enhanced flow through
- Compatible with 6-well plates
- Included adapter for tubes smaller than 50mL



Cat No.	Description	Qty	Price
Cell Strainer, Sterile, Individually Wrapped, w/ One Reducing Adapter			
C7005	5µm Nylon Mesh, Orange PP Frame, Sterile	50	\$120.45
C7010	10µm Nylon Mesh, Teal PP Frame, Sterile	50	\$120.45
C7020	20µm Nylon Mesh, Red PP Frame, Sterile	50	\$120.45
C7030	30µm Nylon Mesh, Purple PP Frame, Sterile	50	\$120.45
CS-4040	40µm Nylon Mesh, Blue PP Frame, Sterile	50	\$120.45
C7050	50µm Nylon Mesh, Grey PP Frame, Sterile	50	\$120.45
CS-4070	70µm Nylon Mesh, White PP Frame, Sterile	50	\$120.45
CS-4100	100µm Nylon Mesh, Yellow PP Frame, Sterile	50	\$120.45



Max RCF

28,000×g

Temperature

Min: -90°C
Max: 121°C


Material

Polypropylene

SureSeal S™ Microtubes

Snap Cap Microcentrifuge Tubes

- Sterile comes standard - E-Beam sterilized
- Self-Standing bag with gusseted bottom
- RNase, DNase, and pyrogen free



Cat No.	Description	Qty	Price
207-GMT*	1.5mL Microtube, Attached Snap Cap, Natural Color, Self-Standing Bag	500	\$18.80
C2000-50	1.5mL Microtube, Attached Snap Cap, Natural Color, 10 Bags of 50	500	\$22.39
C2000-50AST	1.5mL Microtube, Attached Snap Cap, Assorted Colors, 10 Single-Color Bags of 50	500	\$23.51
207-GMTA	1.5mL Microtube, Attached Snap Cap, Assorted Colors, Self-Standing Bag	500	\$21.40
208-GMT	2.0mL Microtube, Attached Snap Cap, Natural Color, Self-Standing Bag	500	\$21.08
202-GMT	0.5mL Microtube, Attached Snap Cap, Natural Color, Self-Standing Bag	500	\$21.08
202-GMTA	0.5mL Microtube, Attached Snap Cap, Assorted Colors, Self-Standing Bag	500	\$22.41

*Also available in individual colors: Blue, Green, Red, Purple, Yellow

SpectraTubes™

Color Cap Centrifuge Tubes

- Color coded for easy sample identification
- Sterile; Available racked or bagged
- 15mL rated to 12,000×g
- 50mL rated to 15,000×g
- RNase, DNase, pyrogen and cytotoxin free



Cat No.	Description	Qty	Price
C2715	15mL Tubes, Sterile, 5 Colors per Foam Rack, 10 racks of 50	500	\$225.51
C2716	15mL Tubes, Sterile, Individual Color Bags, 20 Bags of 25 Tubes, 4 Bags per Color	500	\$205.50
C2750	50mL Tubes, Sterile, 5 Colors per Foam Rack, 20 racks of 25	500	\$238.85
C2751	50mL Tubes, Sterile, Individual Color Bags, 20 Bags of 25 Tubes, 4 Bags per Color	500	\$217.50

ClearSeal™

Screw Cap Microtubes

- Leak-proof O-Ring seal
- Ultra-clear polypropylene
- Resealable sterile bags of 50
- RNase, DNase, and Pyrogen free
- White printed graduations
- Conical or standing bottoms

Temperature
Min: -86°C
Max: 121°CMaterial
PolypropyleneMax RCF
22,000×g

BUY 3 GET 1
FREE!


0.5 mL
Conical0.5 mL
Self-Standing1.5 mL
Conical1.5 mL
Self-Standing2.0 mL
Conical2.0 mL
Self-Standing2.0 mL
Self-Standing
Loop Cap

Cat No.	Description	Qty	Price
C3205-C	0.5mL MicroTube, Sterile, Conical Bottom, O-Ring Screw Cap,	1000	\$271.38
C3205-S	0.5mL MicroTube, Sterile, Self-Standing, O-Ring Screw Cap,	1000	\$271.38
C3215-CG	1.5mL MicroTube, Sterile, Conical Bottom, O-Ring Screw Cap,	1000	\$271.38
C3215-SG	1.5mL MicroTube, Sterile, Self-Standing, O-Ring Screw Cap, Graduated	1000	\$271.38
C3220-CG	2.0mL MicroTube, Sterile, Conical Bottom, O-Ring Screw Cap, Graduated	1000	\$271.38
C3220-SG	2.0mL MicroTube, Sterile, Self-Standing, O-Ring Screw Cap, Graduated	1000	\$271.38
C3220-SGL	2.0mL MicroTube, Sterile, Self-Standing, O-Ring Loop Screw Cap, Graduated	1000	\$283.21

Centrifuge Tubes

5mL to 500mL Capacities

- Vivid graduations
- Perfect size for any application
- RNase, DNase, & pyrogen free

Temperature
Min: -86°C
Max: 121°CTube: Polypropylene
Cap: Polyethylene

5mL

10mL

15mL

25mL

30mL

50mL

50mL
Self-Standing

175mL

225mL

250mL

500mL

Cat No.	Description	Qty	Price	Cat No.	Description	Qty	Price	Cat No.	Description	Qty	Price
5mL MacroTubes®				25mL Centrifuge Tubes				Large Volume Centrifuge Tubes			
T2530	Separate Caps, Non-Sterile, Bulk Bagged	500	\$109.40	C2625-B	Attached Cap, Sterile, Bags of 25	200	\$111.83	C0175	175mL, Sterile PP, Conical Bottom, Bags of 6	120	\$256.39
T2535	Attached Caps, Sterile, Bags of 100	500	\$116.04	C2625-R	Attached Cap, Sterile, Racks of 25	200	\$115.65	C0225	225mL, Sterile, PP, Conical Bottom, Bags of 5	90	\$215.49
T2540	Attached Caps, Sterile, Racks of 50	500	\$124.16	30mL Centrifuge Tubes				C1250	250mL, Sterile, PP, Conical Bottom, Bags of 5	70	\$258.59
T2540-OB	Attached Caps, Opaque Black, Sterile, Racks of 50	500	\$143.11	C2630	Screw Cap, Self-Standing, Sterile, Bags of 25	500	\$186.81	C1500	500mL, Sterile, PP, Conical Bottom, Bags of 5	40	\$264.12
T2545	Attached Caps, Vivid Graduations, Sterile, Racks of 50	500	\$131.50	50mL Centrifuge Tubes				Centrifuge Tube Racks			
10mL Centrifuge Tubes				C2602	Attached Cap, Sterile, Racks of 25	500	\$220.47	R1030	4×50/25mL, 12×15/10/5mL, 32×1.5/0.5mL Tubes	5	\$55.70
C2410	Attached Cap, Sterile, Bags of 50	1000	\$258.59	C2603	Attached Cap, Sterile, Bags of 25	500	\$198.25	T2592	Decanting, 50×5mL, Silicone Grip Wells	4	\$52.38
15mL Centrifuge Tubes				C2603-B	Separate Screw Cap, Non-Sterile, Bulk Bag	500	\$181.61	R1070	30 Slots for 10/15mL, 20 Slots for 25/50mL	5	\$69.77
C2600	Attached Cap, Sterile, Racks of 50	500	\$190.63	C2608	Attached Cap, Sterile, Eco-Plastic Racks of 25	500	\$220.47	R1900	6 Slots for 175, 225, or 250mL Tubes	5	\$76.25
C2601	Attached Cap, Sterile, Bags of 50	500	\$179.03	C2650-W	Attached Cap, Sterile, Individually Wrapped	300	\$158.10	C1500-RK	4 Slots for 500mL Tubes	1	\$95.16
C2607	Attached Cap, Sterile, Eco-Plastic Racks of 25	500	\$190.63	C2605	Attached Cap, Sterile, Self-Standing, Bags of 25	500	\$227.48				
C2615-W	Attached Cap, Sterile, Individually wrapped	300	\$129.82								

Choose one
FREE

With purchase of item below*

(Free with items: BSH6000, C2417, D1036, V0020, H2606, H1001-M)



Value: \$471

Smart phone
imaging enclosure



Value: \$229

Fixed Speed
Vortexer: 2800 rpm



Value: \$441

Tube Rotator
w/ tube holders

Your choice of:
*Promo Code LM-RE2025

SmartDoc™
Item: E5001-SD

Vornado™
Item: BV101-B

RotoMini™
Item: R2020

ISO BLOCK™

Two Independent Block Chambers

- Two isolated chambers with independent temperature control
- Digital control up to 105°C
- Hinged lid improves uniformity
- Quick-Flip universal block accepts 24 x 1.5mL, 48 x 0.2mL, 14 x 0.5mL



Additional blocks available for all common tubes & vials

BSH6000	isoBlock™ Digital Dry Bath, two independent chambers	\$987
BSWCMB	Quick-Flip™ Block, 24 x 1.5mL, 48 x 0.2mL, 14 x 0.5mL	\$212
BSW02	Block, 48 x 0.2mL tubes or 6 PCR strips	\$158

MC-24™ Touch

High Speed Microcentrifuge

- Full color touch screen
- Up to 16,800 xg
- mySpin™ 1 button program recall
- COMBI-Rotor™ accepts tubes & strips
- Optional rotors for 5mL tubes & 2mL cryovials



Included COMBI-rotor:
Holds 24 x 1.5/2.0mL AND
2 x PCR Strips in one rotor!



C2417	MC-24 Touch Microcentrifuge w/ 24 place COMBI-Rotor	\$2,693
C2417-R5	Optional 8 x 5mL rotor for MC-24 Touch	\$1,150
C2417-R5-ADP2	Cryovial adapter for C2417-R5 rotor, 8/pk	\$143



BEADBUG™ 6

MICROTUBE HOMOGENIZER

- Fast, efficient bead homogenization
- Process up to six samples simultaneously
- Sealed tubes - no cross contamination
- Full range bead materials and sizes available to process a wide variety of sample types

D1036	BeadBug 6 Microtube Homogenizer	\$2,983
D1036-A5	5mL adapters for BeadBug 6, set of 2	\$101
D1032-01	Triple-Pure Zirconium, Ø 0.1mm, 2mL	\$132

Additional bead sizes/materials available. See website for details

ASPIRE

LABORATORY ASPIRATOR

- Built-in vacuum pump
- Safe & precise aspiration
- Compact design
- Adjustable vacuum pressure
- Vacuum gauge



V0020	ASPIRE Lab Aspirator with internal pump, 2L PC bottle, tubing, handheld controller & single channel adapters	\$1,108
-------	--	---------

Additional adapters available. See website for details

SURETEMP™

Digital Ovens

- Temperatures up to 200°C
- Large viewing window for monitoring
- Stainless steel chamber
- Fiberglass insulation
- In-lab calibration feature
- Three sizes: 30L, 50L & 150L



H2606-30	SureTherm™ Digital Oven, 30L	\$1,295
H2606-50	SureTherm™ Digital Oven, 50L	\$1,595
H2606-150	SureTherm™ Digital Oven, 150L	\$2,895

INCU-SHAKER™

Mini

- MAGiC™ Clamp accessories for instant-magnetic clamp exchange
- NEW! Adjustable clamp for 125, 250 & 500mL flasks
- Full color, touch screen display
- Compact design with an impressively small footprint of just 11.5 x 9.5 inches



Speed range:	30-300 rpm
Flask capacity:	2 x 1L (H1000-MR-1000) 5 x 500mL (H1000-MR-500) 8 x 250mL (H1000-MR-250) 12 x 125mL (H1000-MR-125)



H1000-MR-CMB
Adjustable clamp for
125, 250, and 500mL flasks

H1001-M	Incu-Shaker Mini, includes standard rubber mat platform	\$3,635
H1000-MR	MAGiC™ Clamp universal platform	\$314
H1000-MR-CMB	MAGiC™ Clamp adjustable clamp For 125mL, 250mL, 500mL flasks	\$56.38