

Premium Fetal Bovine Serum



- Triple 0.1 μm Sterile Filtered
- Free from virus and Mycoplasma
- Endotoxin < 1 EU/ml
- Hemoglobin level < 15 mg/dL
- Heat inactive available

Catalog#	Description	Size	List Price
P711S50	USDA Approved Origin	50ml	\$66.00
P711S500	USDA Approved Origin	500ml	\$520.00
P711U50	US Origin	50ml	\$118.80
P711U500	US Origin	500ml	\$840.00

- Stored at -20°C
- Available in 50ml and 500ml
- No freezing and thawing

Premium Fetal Bovine Serum provides essential nutrients and growth factors to support proliferation, primary cultures, and assays.