



Gobble up these amazing BLACK FRIDAY DEALS!

BLACK FRIDAY SALE

ALL MONTH LONG

BENCHMIXER™ MP

Vortex Mixer for Tubes and Plates

NEW

- * Exceptionally versatile, supplied with:
 - Microplate head (fits a variety of plates)
 - Adapters for 1.5/2.0ml, 15ml and 50ml tubes
 - Cup head to convert to standard vortexer
- * Digital control
 - 200-3,000rpm, up to 999 minutes or continuous

BV1007 BenchMixer MP for microplates and tubes with adapters for 30x1.5/2ml, 12x15ml, 4x50ml tubes and cup head ~~\$499~~ **\$349**

INTRODUCTORY PRICE - SAVE \$150!



Supplied with unit



MC-24R™

REFRIGERATED MICROCENTRIFUGE

- * Touch control and digital display
- * Powerful refrigeration, as low as -20°C
- * High speed - 16,800 xg
- * Quiet operation
- * Supplied complete with COMBI-Rotor™
 - Accepts 24x1.5/2.0ml tubes AND 2xPCR strips



C2417-R MC-24R Refrigerated Microcentrifuge ~~\$5695~~ **\$4195**
 Optional rotors:
 C2417-2420H 24x1.5/2.0ml hermetically sealed rotor \$841
 C2417-3202 4xPCR strip rotor (32x0.2ml) \$795



COMBI-Rotor

SPECIAL BLACK FRIDAY PRICE - SAVE \$1500!

Call Us Today!
BLACK FRIDAY
DEALS END
12/6/2024



TC 9639 Gradient Thermal Cycler

- * Intuitive touch screen control
- * Program wizard - quick, easy and foolproof
- * Multifformat block for 96x0.2ml tubes, PCR strips, 39x0.5ml tubes or a 96 well plate
- * Advanced programming features such as adjustable ramp rate and time/temperature increments/decrements



**SALE!
FREE
PLATEFUGE
WITH
PURCHASE!**

T5000-96 TC 9639 Gradient Thermal Cycler ~~\$5995~~ **\$4995**
 C2000 PlateFuge™ Microplate Microcentrifuge ~~\$894~~ **FREE!**

FREE PLATEFUGE - \$894 VALUE PLUS \$1000 DISCOUNT!





Black Friday BOGO Deals



SMART SAFE Nucleic Acid Stain GLOW

- * Excellent replacement for Ethidium Bromide
 - * Excitation by UV or blue light
 - * Non-mutagenic, environmentally safe
 - * Available as a pre stain (PS) or loading dye (LD)
- Buy one, get one FREE! Mix or Match!
Stock up and SAVE!**

E4500-PS SmartGlow Pre stain, 1ml (20,000X) \$124
E4500-LD SmartGlow Loading Dye with stain, 1ml (6X) \$124



BUY ONE AT \$124 - GET ONE FREE!



BEADBUG MICROTUBE HOMOGENIZER

- * Powerful mixing action for lysis, grinding and homogenization
- * Eliminates cross contamination
- * Personal sized - 3 tube capacity
- * Cold room compatible
- * Wide variety of grinding media available



D1030 BeadBug 3, 3x2ml tube capacity \$998/2 units
D1032-30 TriplePure Zirconium Beads, 3mm, pk 50 tubes \$182
 Additional bead materials/sizes available - see special offer below!

Perfect for small labs and student labs!

BUY ONE AT \$998 - GET ONE FREE!



Benchmark Homogenizer Bead Tubes

- * Convenient, prefilled, 2ml bead tubes
- * Compatible with most brands of bead beaters
- * Variety of bead sizes and materials to meet a range of sample needs
- * TriplePure™ line certified nuclease, protease and DNA free

D1032-30	TriplePure Bead Tubes, Prefilled, 3.0mm zirconium, pk 50 tubes	\$182
D1032-15	TriplePure Bead Tubes, Prefilled, 1.5mm zirconium, pk 50 tubes	\$182
D1032-10	TriplePure Bead Tubes, Prefilled, 1.0mm zirconium, pk 50 tubes	\$182
D1032-05	TriplePure Bead Tubes, Prefilled, 0.5mm zirconium, pk 50 tubes	\$182
D1032-01	TriplePure Bead Tubes, Prefilled, 0.1mm zirconium, pk 50 tubes	\$182
D1032-SK	TriplePure Starter Kit, 10 each of 0.1, 0.5, 1.0, 1.5 and 3.0mm	\$182
D1031-10	Prefilled Bead Tubes, 1.0mm glass, pk of 50 tubes	\$131
D1031-05	Prefilled Bead Tubes, 0.5mm glass, pk of 50 tubes	\$131
D1031-01	Prefilled Bead Tubes, 0.1mm glass, pk of 50 tubes	\$131
D1033-28	Prefilled Bead Tubes, 2.8mm steel, pk of 50 tubes	\$182
D1033-30G	Prefilled Bead Tubes, Garnet Shards, 1x6mm zirconium, pk/50	\$182



BUY 2 PACKS, GET ONE FREE - MIX OR MATCH!

BOGO

More BOGOs and Deals



GET ONE FREE



EZ Pack Agarose Tablets

- * 0.5g tablets eliminate weighing
- * Fast dissolving - just 2 minutes
- * Popular Benchmark Agarose LE
- * Convenient blister pack for safe, easy dispensing

AGLE-1TAB EZ Pack Agarose, 200 tablets (100g) \$177/2 pks



BUY ONE FOR \$177, GET ONE FREE!

Spin microtubes and PCR strips without changing the rotor!



myFUGE™ 12

- * COMBI-Rotor™ for tubes and strips
 - 12x1.5/2.0ml tubes
 - 32x0.2ml PCR tubes
 - 4xPCR strips (8x0.2ml)
- * Twice the capacity of most mini centrifuges
- * Near silent operation
- * Starts and stops with closing and opening the lid

C1012 myFuge 12 Mini Centrifuge with COMBI-rotor
C1008-A5-6 Adapters for 0.5ml tubes, pk of 6

\$457/2 units
\$26

BUY ONE, GET ONE FREE!



COMBI-Rotor

BUY ONE AT \$457, GET ONE FREE!

BactiZapper™

- * Replaces open flame burners for sterilization of loops, needles and tubes
- * Adjustable tilt for ergonomic placement
- * No splatter
- * Optional slide dryer and loop holder



B1001 BactiZapper Tilt Microsterilizer \$576/2 units
 B1000-01 Loop holder attachment \$71
 B1000-03 Slider dryer attachment for 3 slides \$141

BUY ONE BACTIZAPPER - GET ONE FREE!

Roto-Mini™ Plus

- * Variable speed rotation, 5-70rpm
- * Advanced shaking/mixing modes: Mix (Oscillating), Pause, Rock
- * Includes holders for common tubes:
 - 11 x 1.5/0.2ml (pack of 2)
 - 5 x 10-15ml (pack of 2)
 - 5 x 5-7ml (pack of 2)
 - 6 x 50ml (pack of 2)



R2024 Roto-Mini Plus Variable Speed Rotator \$530/2 units
 R2020-RP Mini rocking platform w/rubber mat, 4.5x3" \$114
 R2020-PCR Holder, 96x0.2ml, 12xPCR strips, 1xPCR plate \$88

BUY ONE, GET ONE FREE - \$530 VALUE!

Check out these UNBELIEVABLE DEALS!

ACCURIS Compact Balances



- * Ideal for basic weighing needs
- * Multiple modes: weigh, count, percent
- * LCD display, 0.01g readability
- * Gold-plated, ceramic capacitors
- * Easy, in-lab calibration

Buy one, get one FREE! Mix or match sizes!*

W3300-120	Accuris Compact Balance, 120g capacity	\$327
W3300-300	Accuris Compact Balance, 300g capacity	\$381
W3300-500	Accuris Compact Balance, 500g capacity	\$381

*Least expensive item is free

**BOGO
FREE!**



BOGO - MIX or MATCH - SAVE UP TO \$381

**FREE
myFUGE MINI!**

UV CLAVE Ultraviolet Chamber



- * Improves quality of PCR research by reducing potential contaminants
- * Multiple lamps and reflective chamber walls provide true 360° irradiation;
- * Proprietary UVC transparent shelf, eliminates potential for missed coverage (blind spots)
- * Inactivates DNA, RNA, viruses & bacteria on PCR vessels
- * Adjustable cycle time

B1450	UVClave UltraViolet Chamber (research use only)	\$573
C1008	myFuge Mini Centrifuge with 1.5/2ml tube rotor, 0.2ml strip tube rotor and 0.5ml tube adapters	\$346 FREE!



FREE myFUGE MINI WITH PURCHASE - \$346 VALUE!

PLATEFUGE™ MICROPLATE MICROCENTRIFUGE

- * Quickly spins down droplets in PCR plates
- * Unique, patented swing-out rotor design prevents spillage
- * Half the size of traditional centrifuges
- * Accepts all popular types of PCR plates (skirted, semi, and non) and most microtiter plates



C2000	PlateFuge Microplate Microcentrifuge	\$894/2 units
C2000-A02	Adapter for 0.2ml tubes and strips	\$67
C2000-DWMP	Optional rotor for deepwell plates	\$231



BUY ONE, GET ONE - \$894 VALUE!

Free offers require purchase of designated items at listed price. Offers expire 12/06/2024. Not valid on previous purchases. Special price applies only to item specified. Other capacities or items cannot be substituted. Use promo code TURKEY24.



LabMart® Distribution Companies

JH Berge
(908) 561-1234
jhberge.com
info@labmart.com

LPS
(800) 388-0166
lpsinc.com
info@lpsinc.com

Dot Scientific
(800) 878-1785
dotscientific.com
orders@dotscientific.com