


Solution Reservoirs



Spouts on all four corners for easy pouring in any hands



U-Bottom Micro-trough™ minimizes reagent consumption



Standard 25mL and 50mL supplied with protective lids



5mL - ER-005



10mL - ER-010



25mL - ER-0258



LabMart Brand Solution Reservoirs

- Specially designed for 8-Channel Pipettes
- Conserves valuable reagents
- 75mm inner chamber width




LabMart Brand Solution Reservoirs

	Cat. No.	Description	Previous LPS Part #	Qty
5 mL	ER-005NS	Reservoirs, 5mL, for 8-Channels, non-sterile, bulk pack	(M170500)	300
	ER-005S-1	Reservoirs, 5mL, for 8-Channels, sterile, individually wrapped	(M170501)	100
	ER-005S-5	Reservoirs, 5mL, for 8-Channels, sterile, 5 / sterile bag	(M170505)	200
10 mL	ER-010NS	Reservoirs, 10mL, for 8-Channels, non-sterile, bulk pack	(M180100)	300
	ER-010S-1	Reservoirs, 10mL, for 8-Channels, sterile, individually wrapped	(M180101)	100
	ER-010S-5	Reservoirs, 10mL, for 8-Channels, sterile, 5 / sterile bag	(M180105)	200
25 mL (reduced)	ER-0258NS	Reservoirs, 25mL, for 8-Channels, non-sterile, bulk pack		100
	ER-0258S-1	Reservoirs, 25mL, for 8-Channels, sterile, individually wrapped		100
	ER-0258S-5	Reservoirs, 25mL, for 8-Channels, sterile, 5 / sterile bag		200

LabMart Brand Standard Width Solution Reservoirs

- For 8, 12 Channel Pipettes
- Standard Width - 110mm
- 25mL and standard 50mL supplied with lids



25mL w/ Lid - ER-025



50mL w/ Lid - ER-050



50mL Divided - ER-051



100mL - ER-0100

LabMart Brand Standard Width Solution Reservoirs

	Cat. No.	Description	Previous LPS Part #	Qty
25 mL (standard)	ER-025NS	Reservoirs, 25mL non-sterile, w/ 5 Lids, bulk pack	(M180125)	100
	ER-025S-1	Reservoirs, 25mL sterile, w/ 5 Lids, Individually wrapped	(M180126)	100
	ER-025S-5	Reservoirs, 25mL sterile w/ 5 Lids, 5 / sterile bag	(M180127)	200
	ER-025LID	Lids for 25mL solution reservoirs, sterile, 5 / bag	(M180128)	200
50 mL	ER-050NS	Reservoirs, 50mL non-sterile, bulk pack	(M180150)	100
	ER-050S-1	Reservoirs, 50mL sterile, Individually wrapped	(M180151)	100
	ER-050S-5	Reservoirs, 50mL sterile, 5 / sterile bag	(M180155)	200
	ER-050LID	Lids for 50mL solution reservoirs, sterile, 5 / bag		200
50 mL (divided)	ER-051NS	Reservoirs, 50mL, divided, non-sterile, bulk pack	(M180510)	100
	ER-051S-1	Reservoirs, 50mL, divided, sterile, individually wrapped	(M180511)	100
	ER-051S-5	Reservoirs, 50mL, divided, sterile, 5 / sterile bag	(M180515)	200
100 mL	ER-100NS	Reservoirs, 100mL non-sterile, bulk pack		100
	ER-100S-1	Reservoirs, 100mL sterile, Individually wrapped		100
	ER-100S-5	Reservoirs, 100mL sterile, 5 / sterile bag		200



LabMart® Distribution Companies

JH Berge
(908) 561-1234
jhberge.com
info@labmart.com

LPS
(800) 388-0166
lpsinc.com
info@lpsinc.com

Dot Scientific
(800) 878-1789
dotscientific.com
orders@dotscientific.com